custom walker adjustments

With your arm extended down, handle should be at wrist height. With your hand placed on the hand-grip, your elbow should be flexed 20°-30°.

CAUTION:

Be aware of the following:

- Electrical cords in way
- Liquids on floor
- Throw rugs
- Toys or debris
- Furniture
- Doors and door steps
- Gaps in floor

WARNING:

- Use wheels on front legs only or on all four legs with appropriate rear wheel brakes.
- Do not push one side of the walker ahead and then other side without lifting the walker first.
- Use Drive accessories only. Use of other manufacturers' accessories can result in improper fit and risk fall and/or injury.
- 4. Do not use on stairs or escalators.

lifetime limited warranty

Your Drive branded product is warrantied to be free of defects in materials and workmanship for the lifetime of the product for the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



© 2014 Medical Depot, Inc. All rights reserved.

Drive is a trademark of Medical Depot, Inc.

Port Washington N.Y. 11050 USA Made in China

steel folding walker two button with 5" wheels junior

item # 10247-4



 ϵ

EU Authorized Representativ Drive Medical LTD Ainley's Industrial Estate Elland, West Yorkshire, United Kingdom HX5 9JP

www.drivemedical.com

walker instructions



REF. 10247-4 Shown

Weight Capacity: 350 lbs (158 Kg)



wheelchair to walker transfer

Making certain wheels on the wheelchair are locked lift your body by pushing down on the armrests and carefully transfer one hand at time to the walker.





wheelchair to walker transfer (with lower extremity injury to right leg or foot)

Making certain wheels on the wheelchair are locked, lift your body by pushing down on the armrests and carefully transferring first the left hand, then the right hand to the walker. (If injury to left leg or foot transfer the right hand first, followed by left.)





three-point gait

(for walking with walker)

Step 1

Align the middle of your foot with the back legs of the walker.



Step 2

Slide the walker forward and place it comfortably in front of you.



Step 3

Carefully, walk forward one step at a time.

