



Select store Location – AADL Number  
 Calgary - AADL 00000 78000 003  
 Edmonton- AADL 00000 78000 001  
 Grande Prairie - AADL 0000 665538 001  
 Red Deer - AADL 0000 637100 001

# STELLAR HD

AADL Cat. # W321

**Order Form - AADL**  
 Prices in Canadian Dollars. October 1, 2019

## Customer and Order Information

Purchase Order *	Tag	v.0301
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <span style="float: right;">*required</span>

## Frame Type

<input checked="" type="checkbox"/> 2101000	<b>Stellar HD 45° Bariatric Manual Tilt Wheelchair</b> Extended delivery lead times may apply.	W321	\$ 3,695
<input checked="" type="checkbox"/> 2101201	<b>Reinforced Steel Frame Package</b> Heavy-gauge steel reinforcement of frame and other critical components.		\$ 1,295
<b>i</b> Consider the <b>Stellar Impact</b> for highly agitated users up to 250 lbs. See separate order form.			

## Frame Options

<b>Heavy Duty Base and Seat Frame</b>			
<input type="checkbox"/> 2102045	450 lbs Weight Capacity. 45° of Tilt.		
<input type="checkbox"/> 2102060	600 lbs Weight Capacity. 25° of Tilt.		\$ 300
<b>User Weight Range</b> Required for selection of Gas Struts (Force measured in Newton - N)			
<input type="checkbox"/> 2102212	Up to 250 lbs (2 x 1200 N)	<input type="checkbox"/> 2102308	250 lbs to 350 lbs (3 x 800 N)
<input type="checkbox"/> 2102310	350 lbs to 425 lbs (3 x 1000 N)	<input type="checkbox"/> 2102312	425 lbs to 600 lbs (3 x 1200 N)

## Seat Width Overall chair width is +11"

<b>Standard Seat Widths</b>			
<input type="checkbox"/> 2103022	22" Seat Width	<input type="checkbox"/> 2103024	24" Seat Width
<input type="checkbox"/> 2103028	28" Seat Width	<input type="checkbox"/> 2103030	30" Seat Width
<input type="checkbox"/> 2103026	26" Seat Width	<input type="checkbox"/> 2103032	32" Seat Width
			\$ 350
			\$ 750

## Seat Depth

<b>Standard Seat Depths (Adjustable)</b>			
<input type="checkbox"/> 2104018	18" Seat Depth	<input type="checkbox"/> 2104019	19" Seat Depth
<input type="checkbox"/> 2104021	21" Seat Depth	<input type="checkbox"/> 2104022	22" Seat Depth
<input type="checkbox"/> 2104020	20" Seat Depth		
			\$ 250

## Frame Color

<b>PDG Premium Colors</b>			
<input type="checkbox"/> 2105000	Tracer Red	<input type="checkbox"/> 2105001	Gloss Black
<input type="checkbox"/> 2105003	Black Shimmer	<input type="checkbox"/> 2105004	Orange Fusion
<input type="checkbox"/> 2105006	Blue Shimmer	<input type="checkbox"/> 2105007	Gloss White
<input type="checkbox"/> 2105009	Geyser Green	<input type="checkbox"/> 2105010	Lollypop Blue
<input type="checkbox"/> 2105012	Bubble Gum Pink	<input type="checkbox"/> 2105013	Berry Purple
<input type="checkbox"/> 2105002	Silver	<input type="checkbox"/> 2105005	Karl's Imagination
<input type="checkbox"/> 2105008	Dark Titanium	<input type="checkbox"/> 2105011	Lollypop Red
<input type="checkbox"/> 2105014	Safety Yellow		
<input type="checkbox"/> 2105200	<b>Custom Color</b> Contact PDG for specifications.		\$ 300

## Arm Type

<input type="checkbox"/> 2106000	<b>No Armrest</b>	
<input type="checkbox"/> 2106214	<b>Dual Post Adjustable Height (10"-14")</b>	\$ 275

## Armrests

### Regular Armpads

- 2107100 Desk Length - 2" x 10"       2107101 Full Length - 2" x 14"

### Wide Armpads

May interfere with wheels in some tilt configurations.

\$ 200

- 2107110 Desk Length - 3" x 10"       2107111 Full Length - 3" x 12"

### Gel Armpads

May interfere with wheels in some tilt configurations.

\$ 135

- 2107120 Full Length - 2" x 12"       2107121 Full Length - 3.5" x 12"

### Contour Arm Trough with Hand Pad

May interfere with wheels in some tilt configurations.

\$ 125

- 2107130 Right (4.5" x 18")       2107131 Left (4.5" x 18")

### Clothing Guards

\$ 165

## Seat Pan

- 2108001 **Solid Aluminum Seat Pan**

\$ 155

## Cushion

## Seat Height

### Seat to Floor Height

- 2110016 Set at 16" STF       2110017 Set at 17" STF       2110018 Set at 18" STF

## Casters

- 2111050 **5" Caster With Heavy Duty Steel Fork**

## Rear Wheels

- 2112224 **24" Heavy Duty Spoke**

## Axle Type

- 2113010 **Fixed Rear Axles**

## Rear Tires

- 2114000 **Urethane**

## Hand Rims

- 2115000 **No Hand Rims**

- 2115010 **Standard**

- 2115020 **Plastic Coated**

\$ 250

## Wheel Locks

- 2116010 **Push to Lock**

- 2116020 **Pull to Lock**

- 2116100 **Wheel Lock Extensions**

\$ 60

## Anti-Tippers

- 2117010 **Rear Anti-Tippers** Required Safety Option.

\$ 60

## Back Post Style

<input checked="" type="checkbox"/>	2118010	<b>Push Bar Back</b>	
<b>Back Post Height</b> Seat Pan to Handle			
<input type="checkbox"/>	2118425	25" Height	<input type="checkbox"/> 2118431 31" Height

## Back Options

<input checked="" type="checkbox"/>	2119090	<b>Seat to Back Angle set at 90°</b>	
<input type="checkbox"/>	2119240	<b>Dynamic Rocker Back Assembly</b>	\$ 2,995
Not compatible with 30° Recline. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.			

## Back Upholstery

<input type="checkbox"/>	2120000	<b>No Back Upholstery</b>	
<b>Adjustable Tension</b> Back height measured from Seat Pan to top of Upholstery.			
<input type="checkbox"/>	2120321	Set at 21"	<input type="checkbox"/> 2120322 Set at 22" <input type="checkbox"/> 2120323 Set at 23" <span style="float: right;">\$ 295</span>

## Positioning Backrest

## Tilt/Recline Options

<b>Attendant Tilt only</b>			
<input type="checkbox"/>	2123011	Trigger on Right Back Post	<input type="checkbox"/> 2123012 Trigger on Left Back Post
<b>Tilt Limit</b> Maximum Tilt Limit for 600 lbs Weight Capacity is 25°.			
<input type="checkbox"/>	2123220	Set at 20°	<input type="checkbox"/> 2123225 Set at 25° <input type="checkbox"/> 2123230 Set at 30°
<input type="checkbox"/>	2123235	Set at 35°	<input type="checkbox"/> 2123240 Set at 40° <input type="checkbox"/> 2123245 Set at 45°
<input type="checkbox"/>	2123310	<b>Reclining Backrest (Manual)</b>	30° maximum recline. Not compatible with Tension Adjustable Upholstery. <span style="float: right;">\$ 895</span>

## Headrest

<input type="checkbox"/>	2124010	<b>Headrest Mounting Plate</b>	Not required if mounting to a solid Back. <span style="float: right;">\$ 95</span>
--------------------------	---------	--------------------------------	--

## Front Rigging

<input type="checkbox"/>	2125001	<b>Sector Blocks only</b>	
<input type="checkbox"/>	2125130	<b>Swing Away 60°</b>	
<input type="checkbox"/>	2125200	<b>Swing Away Elevating Legrest</b>	Includes Calf Pads. <span style="float: right;">\$ 375</span>

## Foot Plates

<input checked="" type="checkbox"/>	2126020	<b>Aluminum Angle Adjustable</b>	<span style="float: right;">\$ 175</span>
-------------------------------------	---------	----------------------------------	---

## Footrest/Legrest Accessories

<input type="checkbox"/>	2127030	<b>Swing Away ELR Gel Knee Pad</b>	<span style="float: right;">\$ 100</span>
<input type="checkbox"/>	2127040	<b>Adjustable Length Heel Loops</b>	<span style="float: right;">\$ 50</span>

## Extension Tubes

<input checked="" type="checkbox"/>	2128010	<b>Regular Extension Tubes (13"-17")</b>	
-------------------------------------	---------	--	--

## Positioning Straps

--	--	--	--

## Miscellaneous Options

<input type="checkbox"/> 2130120	<b>O2 Holder</b>	\$ 395
<input type="checkbox"/> 2130170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.	\$ 10

## Order Notes



Please use the AADL number of the Eco Medical Store you are working with

Calgary  
AADL 00000 78000 003

Edmonton  
AADL 00000 78000 001

Grande Prairie  
AADL 0000 665538 001

Red Deer  
AADL 0000 637100 001