

Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001



Yes, you can.®

Invacare® 9000XT Recliner Series
AADL Customers Only
AADL # W615
Effective October 1, 2021

NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

FRAME & SEAT

CHAIR TYPE

- 9RC 9000XT Recliner \$1,619

FRAME WIDTH & DEPTH

- WD46 14" Wide X 16" Deep No Charge
- WD47 14" Wide X 17" Deep \$145
- WD48 14" Wide X 18" Deep \$145
- WD66 16" Wide X 16" Deep No Charge
- WD67 16" Wide X 17" Deep \$145
- WD68 16" Wide X 18" Deep \$145
- WD86 18" Wide X 16" Deep No Charge
- WD87 18" Wide X 17" Deep \$145
- WD88 18" Wide X 18" Deep \$145
- WD96 19" Wide X 16" Deep No Charge
- WD97 19" Wide X 17" Deep \$145
- WD98 19" Wide X 18" Deep \$145
- WD06 20" Wide X 16" Deep No Charge
- WD07 20" Wide X 17" Deep \$145
- WD08 20" Wide X 18" Deep \$145
- WD26 22" Wide X 16" Deep \$206
- WD27 22" Wide X 17" Deep \$333
- WD28 22" Wide X 18" Deep \$333

FRAME TYPE

- R Recliner \$774

This is an automatic pick in the configurator.

ARM & BACK

ARM TYPE

Fixed Height Arms 9", Adjustable Height Arms 9"-13", removable style.

- 18 Full Length Fixed Height Conventional \$261
- 28 Desk Length Fixed Height Conventional \$261
- 19 Full Length Adj. Height Conventional \$372
- 29 Desk Length Adj. Height Conventional \$372
- 19JY Full length Adj. Height Conventional (8"-10") \$372
- 29JY Desk length Adj. Height Conventional (8"-10") \$372

BACK TYPE

Add additional 11" for the headrest extension.

- 2R Reclining Back 23" high No Charge
- L Reclining Back 23" high w/Lift Package \$130

COLOURS

FRAME FINISH

- 17P Cobalt Blue No Charge
- 61P Electric Red No Charge
- 24P Wet Black No Charge
- 34P Satin (Flat) Black **Standard**
- 71P Silver Vein No Charge

CHAIR FEATURES AND OPTIONS

- ◆ Chrome plated steel frame
- ◆ Choice of frame colours available
- ◆ Conventional adjustable and fixed-height arm styles
- ◆ 24" composite/Pneumatic tire
- ◆ Degree of recline from 90° to 168°
- ◆ 22" high back and 10" headrest extension
- ◆ {removable} with pillow standard
- ◆ Non-swing back removable arm
- ◆ Attendant operated wheel locks
- ◆ Lifetime warranty on frame and crossbraces
- ◆ Standard with rear anti-tippers
- ◆ STF = 19.5" with 24" wheels
- ◆ STF = 18.5" with 22" wheels
- ◆ STF = 17.5" with 20" wheels

WEIGHT LIMITATION 250 LBS.

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

UPHOLSTERY COLOURS

All back upholstery is screwed on.

- U240 Black Nylon Upholstery **Standard**
- U245 Black Sail Cloth Adjust. Tension Seat **No Charge**
- U245 sail cloth adjustable tension seat upholstery available w/U240 black nylon back upholstery only.*
- U550 Midnight Blue Vinyl Upholstery **No Charge**
- U67 Black Vinyl **No Charge**
- U75 Dove Grey Vinyl **No Charge**
- U86 Ruby Vinyl **No Charge**
- U88 Camel Vinyl **No Charge**

CASTERS

CASTER OPTIONS

- 1228 8" X 1 3/4" Semi-Pneu Caster Lt. Grey **\$150**
- AT903 8" X 1-1/4" Composite/Urethane **Standard**
- 1290 8" X 1" Solid Rubber Casters Charcoal **No Charge**
- 1299 8" X 1" Solid Rubber Casters Dark Grey **No Charge**
- 1292 8" X 1-1/4" Pneumatic Casters **No Charge**
- 1492 8" X 1-1/4" Pneu w/Flat Free Insert **\$128**

Seat to Floor Height	Seat-to-Floor Matrix 8" Casters		
	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17.50"	Top	20"	1 Position Only
18.50"	Middle	22"	1 Position Only
19.50"	Bottom	24"	1 Position Only

WHEELS

WHEEL OPTIONS

- 1763BK 20" Composite/Pneumatic Tires **\$231**
- 4763BK 20" Composite/Pneumatic-Flat Free Inserts **\$349**
- U1763BK 20" Composite/Urethane Tires **No Charge**
- 1764BK 22" Composite/Pneumatic Tires **\$231**
- 4764BK 22" Composite/Pneumatic-Flat Free Inserts **\$349**
- U1764BK 22" Composite/Urethane Tires **No Charge**
- 2244BK 24" Composite/Pneumatic Tires **Standard**
- 4244BK 24" Composite/Pneumatic-Flat Free Inserts **\$128**
- U2244BK 24" Composite/Urethane Treaded Tires **No Charge**
- U2222C 24" Composite/Urethane Corded Tires **No Charge**

Not available with AHDRM or 1225-12.

HANDRIM OPTIONS

- AHDRM Aluminum Handrims **Standard**
- Not available with U2222C rear wheels.*
- COM Composite Handrims **No Charge**
- Available with 24" rear wheels only.*
- 1512 Plastic Coated Aluminum Handrims **\$114**
- 1224-8 8 Alum Oblique Proj. Hdrm. Plastic Coated **\$307**
- Available with 24" rear wheels only.*
- 1225-8 8 Alum Vert. Proj. Hdrm., Plastic Coated **\$307**
- Available with 24" rear wheels only.*

WHEEL LOCKS

- Attendant Operated **Standard**
- 1329BK Toggle Wheel Lock Extension (8") **\$63**
- 9696R Push to Lock Wheel Locks (Patient Operated) **\$48**
- Attendant operated wheel locks are standard and are not available with patient operated wheel locks.*
- 1245 Pull to Lock Wheel Locks (Patient Operated) **No Charge**
- Attendant operated wheel locks are standard and are not available with patient operated wheel locks.*

POSITIONING STRAPS

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap **\$84**
- 1515 Auto Style Buckle Seat Positioning Strap **\$84**
- 1321BK Hook & Loop Seat Positioning Strap **\$66**

CHEST POSITIONING STRAPS

- 9515 Auto Style Chest Positioning Strap **\$84**



Select store location - AADL Numbers

- Calgary - AADL 0000 78000 003
- Edmonton - AADL 0000 78000 001
- Grande Prairie - AADL 0000 665538 001

ACCESSORIES & OPTIONS

BACK ACCESSORIES

<input type="checkbox"/>	1160	Carrying Pocket	\$83
<input type="checkbox"/>	3040	Torso Support for Recliner	\$175
<input type="checkbox"/>	STC	Stencil Back Upholstery	\$59
<i>3/4" block letters. Only available with vinyl upholstery.</i>			
<input type="checkbox"/>	3046	Head Support for Recliner	\$163
<input type="checkbox"/>	RB	Dual Back Liner	\$91
<input type="checkbox"/>	L8601	Omit Back Upholstery	Less Net (\$20)

SEAT ACCESSORIES

<input type="checkbox"/>	9120	Solid Seat Insert	\$144
<input type="checkbox"/>	RS	Dual Seat Liner	\$90

ARM ACCESSORIES

<input type="checkbox"/>	NRAS	Non-Removable Arm With Screws	No Charge
<input type="checkbox"/>	1300	Utility Tray	\$268

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	2219	IV Rod, Telescopic (Recliner)	\$163
<input type="checkbox"/>	2226	Front Anti-Tippers	\$158
<input type="checkbox"/>	1100	Tool Kit	No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING

<input type="checkbox"/>	AHL4A	Hemi Smartleg Artic. Legrests 16" wide Alum. Footplates.	\$399
<i>Legrest available in black finish only.</i>			
<input type="checkbox"/>	AHL4C	Hemi Smartleg Artic. Legrests 16" wide Comp. Footplates.	\$399
<i>Legrest available in black finish only.</i>			
<input type="checkbox"/>	AHS4A	Hemi Smartleg Articulating Legrests 14"-15" Wide	\$392
<i>Legrest available in black finish only. Only available on 14" wide.</i>			
<input type="checkbox"/>	T93HA	Hemi Sw Footrests Aluminum Footplates.	\$159
<input type="checkbox"/>	T93HC	Hemi Sw Footrests Composite Footplates.	\$159
<input type="checkbox"/>	T94HA	Hemi Sw Elevating Legrests Aluminum Footplates.	\$309
<i>Legrest available in black finish only.</i>			
<input type="checkbox"/>	T94HC	Hemi Sw Elevating Legrests Composite Footplates.	\$302
<i>Legrest available in black finish only.</i>			
<input type="checkbox"/>	T94HS	Hemi Sw Elevating Legrests Aluminum Footplate 14" wide	\$158
<i>Legrest available in black finish only. Only available on 14" wide.</i>			
<input type="checkbox"/>	T93HS	Hemi Sw Footrests Aluminum Footplates. 14" wide	\$158
<i>Only available on 14" wide.</i>			
<input type="checkbox"/>	T93HAB	Hemi Bolt-Thru First, Aluminum Footplates	\$158
<input type="checkbox"/>	T93HCB	Hemi Bolt-Thru First, Composite Footplates	\$158
<input type="checkbox"/>	T94HAB	Hemi Bolt-Thru Legrest, Aluminum Footplates	\$302
<input type="checkbox"/>	T94HCB	Hemi Bolt-Thru Legrest, Composite Footplates	\$302
<input type="checkbox"/>	T93HSB	Hemi Bolt-Thru First, Aluminum Footplates 14" wide	\$158
<i>Only available on 14" wide.</i>			
<input type="checkbox"/>	T94HSB	Hemi Bolt-Thru Legrest, Aluminum Footplates 14" wide	\$158
<i>Only available on 14" wide.</i>			

FOOTPLATE OPTIONS

<input type="checkbox"/>	1350	Extra Large Aluminum Footplates	\$106
<input type="checkbox"/>	AT5543	Adjustable Angle Flip-up Footplates	\$132

FOOTREST / LEGREST ACCESSORIES

<input type="checkbox"/>	1820	Non-Removable Device	\$41
<i>Not Available with AHL4A or AHL4C elevating legrests.</i>			
<input type="checkbox"/>	1337	Calf Strap	\$63
<input type="checkbox"/>	1600BK	Heel Loops With Ankle Straps	\$84
<input type="checkbox"/>	1348	Legrest Panel	\$93
<i>Not available with Hemi minimum Swingaway footrests.</i>			
<input type="checkbox"/>	1346	Heel Loops	\$63
<i>Standard on footrests.</i>			
<input type="checkbox"/>	L8689	Omit Heel Loops	Less Net (\$10)
<input type="checkbox"/>	1319	Black Non-Padded Calf Pads	Less Net (\$10)
<i>Not available with AHL4A or AHL4C elevating legrests.</i>			

Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001



Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)
 Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
 WWW.INVACARE.CA