



- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001



Freestander - AADL

September 2020



Mark For: Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: PV: PC: _____ Dealer Phone: Fax: _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> E-Mail	Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ <div style="text-align: center; border: 1px solid black; padding: 2px;">ADDITIONAL SHIPPING INFORMATION</div> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: PV: PC: _____ Ship To Phone: Fax: _____
--	---

Freestander

Freestander
Selection required. Includes Chest and Hip Pad, Flexible Lateral, Lap Belt, PU Knee Supports, Footplate, Sandals, Tray Cover and Activity Bowl.

Description	Item Code	Retail
<input type="checkbox"/> Size 1 Freestander - Meadow	LFS/1-MW-EXP	\$3,140.00
<input type="checkbox"/> Size 2 Freestander - Mondeo Blue	LFS/2-MB-EXP	\$3,400.00
<input type="checkbox"/> Size 3 Freestander - Mondeo Blue	LFS/3-MB-EXP	\$3,847.00

Description	Item Code	Retail
<input type="checkbox"/> Size 1 Lap Belt - Meadow	LFS/1/B-MW	\$227.00
<input type="checkbox"/> Size 1 Lap Belt - Mondeo Blue	LFS/1/B-MB	\$232.00
<input type="checkbox"/> Size 2 Lap Belt - Meadow	LFS/2/B-MW	\$321.00
<input type="checkbox"/> Size 2 Lap Belt - Mondeo Blue	LFS/2/B-MB	\$236.00
<input type="checkbox"/> Size 3 Lap Belt - Mondeo Blue	LFS/3/B-MB	\$256.00
<input type="checkbox"/> Size 1 High Back Strap - Meadow	LFS/1/D-MW	\$142.00
<input type="checkbox"/> Size 2 High Back Strap - Mondeo Blue	LFS/2/D-MB	\$151.00
<input type="checkbox"/> Size 3 High Back Strap - Mondeo Blue	LFS/3/D-MB	\$169.00
<input type="checkbox"/> Size 1 Solid Pelvic Plate - Meadow <i>Pelvic Band selection required</i>	LFS/1/L-MW	\$369.00
<input type="checkbox"/> Size 2 Solid Pelvic Plate - Mondeo Blue <i>Pelvic Band selection required</i>	LFS/2/L-MB	\$416.00
<input type="checkbox"/> Size 3 Solid Pelvic Plate - Mondeo Blue <i>Pelvic Band selection required</i>	LFS/3/L-MB	\$466.00
<input type="checkbox"/> Size 1 Pelvic Band - Meadow	LPS/1/B-MW	\$240.00
<input type="checkbox"/> Size 2 Pelvic Band - Mondeo Blue	LPS/2/B-MB	\$263.00
<input type="checkbox"/> Size 3 Pelvic Band - Mondeo Blue	LPS/3/B-MB	\$285.00

Description	Item Code	Retail (Each)
<input type="checkbox"/> Size 1 Upholstered Knee Cups - Meadow	LPS/1/K-MW	\$147.00
<input type="checkbox"/> Size 2 Upholstered Knee Cups - Mondeo Blue	LPS/2/K-MB	\$147.00
<input type="checkbox"/> Size 3 Upholstered Knee Cups - Mondeo Blue	LPS/3/K-MB	\$147.00
<input type="checkbox"/> Size 4 Upholstered Knee Cups - Mondeo Blue	LPS/4/K-MB	\$147.00
Size 1 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/1/F \$75.00
Size 2 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/2/F \$75.00
Size 3 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/3/F \$75.00

Additional Foot Positioning

Sandal selection required with Sandal Raiser

Description	Item Code	Retail (Each)
<input type="checkbox"/> Small Sandals with Straps (Pair)	152-1600	\$222.00
<input type="checkbox"/> Medium Sandals with Straps (Pair)	152-2600	\$236.00
<input type="checkbox"/> Large Sandals with Straps (Pair)	152-3600	\$255.00
<input type="checkbox"/> Size 1 Sandal Raiser (Each)	Quantity: _____ (1-6) LFS/1/J	\$86.00
<input type="checkbox"/> Size 2 Sandal Raiser (Each)	Quantity: _____ (1-6) LFS/2/J	\$86.00
<input type="checkbox"/> Size 3 Sandal Raiser (Each)	Quantity: _____ (1-6) LFS/3/J	\$86.00

Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001

Headrest

	Description	Item Code	Retail
<input type="checkbox"/>	Headrest Interface Frame <i>Required with Contoured Headrest</i>	LIF/1	\$635.00
	Headrest Bracket - Left and Right	LIF/7	\$140.00
<input type="checkbox"/>	Size 1 Contoured Headrest - Meadow	LIF/2/A-MW	\$428.00
<input type="checkbox"/>	Size 1 Contoured Headrest - Mondeo Blue	LIF/2/A-MB	\$479.00
<input type="checkbox"/>	Size 1 Contoured Headrest - Vinyl	LIF/2/A	\$479.00
<input type="checkbox"/>	Size 2 Contoured Headrest - Mondeo Blue	LIF/2/B-MB	\$431.00
<input type="checkbox"/>	Size 2 Contoured Headrest - Vinyl	LIF/2/B	\$479.00
<input type="checkbox"/>	Size 3 Contoured Headrest - Mondeo Blue	LIF/2/C-MB	\$431.00
<input type="checkbox"/>	Size 3 Contoured Headrest - Vinyl	LIF/2/C	\$479.00

Accessories

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Grab/Suspension Rail	LPS/1/C	\$242.00
<input type="checkbox"/>	Size 2 Grab/Suspension Rail	LPS/2/C	\$242.00
<input type="checkbox"/>	Size 3 Grab/Suspension Rail	LPS/3/C	\$242.00

Technical Specifications

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55.1" / 72.8"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Pad Height (Min / Max)	20.7" / 33.5"	26.8" / 43.3"	38.4" / 56.1"
Width Between Chest Supports (Min / Max)	6" / 10"	8" / 12"	10" / 14"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	800 mm / 1151 mm	1100 mm / 1448 mm	1400 mm / 1849 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Pad Height (Min / Max)	470 mm / 680 mm	680 mm / 975 mm	975 mm / 1350 mm
Width Between Chest Supports (Min / Max)	152 mm / 254 mm	203 mm / 305 mm	254 mm / 356 mm
Activity Tray	400 mm x 505 mm	450 mm x 575 mm	500 mm x 625 mm



Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3 • Concord, ON • Canada L4K 4V3