



Vendor Name: _____
 PO #: _____

Contact Name: _____
 Tag Name: _____ Date: _____

Fax to 905.564.5424 or Email to info@mapleleafwheelchair.ca

Swift - Standard Features

SW-A **\$1895.00**

- ◆ Steel Frame Construction**
- ◆ Standard Back Height - 18"**
- ◆ Removable Height Adjustable Armrests**
- ◆ Weight Capacity - 250 lbs.**
- ◆ Removable Swing Away Footrests**
- ◆ Anti-Tippers**

FRAME COLOUR	CHAIR SIZE				SEAT HEIGHT		ARM TYPE		FOOTREST	
FCBK - Black <input type="checkbox"/>	<u>WIDTH</u>		<u>DEPTH</u>		<u>HEMI FRAME</u>		<i>C-Style = Conventional Arm A-Style = Space Saver</i>		F1PIN60	60° Pin Style <input type="checkbox"/>
	SWCW14	14" <input type="checkbox"/>	SWCD16	16" <input type="checkbox"/>	SWSH1450	14.50" <input type="checkbox"/>	Full Length	A3 <input type="checkbox"/>	F1PIN70	70° Pin Style <input type="checkbox"/>
	SWCW15	15" <input type="checkbox"/>	SWCD17	17" <input type="checkbox"/>	SWSH1525	15.25" <input type="checkbox"/>			Adjustable	C3 <input type="checkbox"/>
	SWCW16	16" <input type="checkbox"/>	SWCD18	18" <input type="checkbox"/>	SWSH1550	15.50" <input type="checkbox"/>	Height	A4 <input type="checkbox"/>	F1PIN90	90° Pin Style <input type="checkbox"/>
	SWCW17	17" <input type="checkbox"/>	SWCD20	20" <input type="checkbox"/>	SWSH1625	16.25" <input type="checkbox"/>			Desk Length	C4 <input type="checkbox"/>
	SWCW18	18" <input type="checkbox"/>		\$209	SWSH1650	16.50" <input type="checkbox"/>	Adjustable	C3L <input type="checkbox"/>	*A2	
	SWCW19	19" <input type="checkbox"/>	SWCD22	22" <input type="checkbox"/>	SWSH1700	17.00" <input type="checkbox"/>	Height			C4 <input type="checkbox"/>
	SWCW20	20" <input type="checkbox"/>		\$350	\$350	SWSH1750	17.50" <input type="checkbox"/>	Full Length	C3L <input type="checkbox"/>	
	UPHOLSTERY		SWCW21	21" <input type="checkbox"/>	<u>STANDARD FRAME</u>		Adjustable			Height Low
	*Seat Upholstery Included	SWCW22	22" <input type="checkbox"/>	*HD Package required on all chairs over 22" Please use Swift W272 Order Form		SWSH1775-S	17.75" <input type="checkbox"/>	Requires 20" or 22" Wheels	ALF	
NU - Nylon <input type="checkbox"/>		\$450	SWSH1825			18.25" <input type="checkbox"/>	SWSH1925		19.25" <input type="checkbox"/>	AAF-A
VU - Vinyl <input type="checkbox"/>			SWSH1975	19.75" <input type="checkbox"/>	SWSH1975	19.75" <input type="checkbox"/>	XLF	Extra Large <input type="checkbox"/>		
SBU - Sling Back Upholstery <input type="checkbox"/>		\$40	SWSH2000	20.00" <input type="checkbox"/>	*Select Wheel Combination Below			Omit Heel Loops <input type="checkbox"/>	(Mee Hn-Middle) \$130	
UBL - Black <input type="checkbox"/>										
WEIGHT					BACK OPTION		BRAKE EXT.		BRAKE TYPE	
Client's Weight: _____ Heavy Duty Package required for client's weight over 250 lbs. Please use Swift W272 order form.					AABK <input type="checkbox"/> Angle Adjustable Back <input type="checkbox"/> \$240 <i>Requires C-style Arms!</i>		BRX6 6" Brake <input type="checkbox"/> Extensions		APUL Aluminum <input type="checkbox"/> Pull To Lock	
							BRX9 9" Brake <input type="checkbox"/> Extensions \$60		APUS Aluminum <input type="checkbox"/> Push To Lock	
									SPUL Steel Pull <input type="checkbox"/> To Lock	
									SPUS Steel Push <input type="checkbox"/> To Lock	

REAR WHEEL <input type="checkbox"/> Urethane - Standard <input type="checkbox"/> Spoke - N/C <small>*not available with HD</small>	CASTER <input type="checkbox"/> x 1 1/4" - N/C <input type="checkbox"/> x 1 3/4" - \$80 <input type="checkbox"/> x 2" - \$80* <small>*only 6" & 8" n/c on HD2/3"</small>	SEAT HEIGHT <small>*STF to the rail</small>	SEAT DUMP	FORK	WHEEL & BRAKE OPTIONS		
Hemi Frame					22AT	22" x 1 3/8" Air Tire	\$88.00 <input type="checkbox"/>
20" x 1 3/8"					24AT	24" x 1 3/8" Air Tire	\$88.00 <input type="checkbox"/>
					PCH	Plastic Coated Handrim	\$106.00 <input type="checkbox"/>
22" x 1 3/8"					AH	Attendant Hand Brake *A1	\$240.00 <input type="checkbox"/>
					UNIBRX	Uni Lateral Wheel Lock *A1 <input type="checkbox"/> Right <input type="checkbox"/> Left	\$117.00 <input type="checkbox"/>
24" x 1 3/8"					QRX	Quick Release Axle	\$87.00 <input type="checkbox"/>
					MLCTD	Transit Tie-Downs *A1 **Unoccupied** WC-19 - up to 250 lbs.	\$160.00 <input type="checkbox"/>
Standard Frame							
22" x 1 3/8"					<input type="checkbox"/> 17.75"	1.00"	3-Hole
					<input type="checkbox"/> 19.75"	0.75"	5-Hole
24" x 1 3/8"					<input type="checkbox"/> 19.25"	1.00"	5-Hole
					<input type="checkbox"/> 18.25"	1.00"	3-Hole
24" x 1 3/8"					<input type="checkbox"/> 19.00"	0.75"	3-Hole
					<input type="checkbox"/> 19.75"	0.75"	5-Hole
<input type="checkbox"/> 20.00"	1.00"						
<input type="checkbox"/> 20.75"	1.00"						



ADDITIONAL OPTIONS

ACCESSORIES				MODIFICATIONS			
NPLS-A	Nylon Padded Leg Strap	\$60.00	<input type="checkbox"/>	CM1	Custom Back Height/Push Handle___ (*6)	\$150.00	<input type="checkbox"/>
O2H-A	Oxygen Tank Holder - Generic	\$130.00	<input type="checkbox"/>	CM2	Special Back Height Upholstery___	\$40.00	<input type="checkbox"/>
IVA	I.V. Pole - Adjustable *A1	\$130.00	<input type="checkbox"/>	CM3	Reinforced Back Post	\$60.00	<input type="checkbox"/>
IVNA	I.V. Pole - Non-Adjustable *A1	\$110.00	<input type="checkbox"/>	CM6-A	1-Piece Flip Up Footboard	\$80.00	<input type="checkbox"/>
ASB-A	Auto Seat Belt - Generic	\$22.00	<input type="checkbox"/>	CM14	Reinforced F/R Bracket	\$60.00	<input type="checkbox"/>
VSB-A	Velcro Seat Belt - Generic	\$15.75	<input type="checkbox"/>	CM15	Reinforced Footrest Pins	\$60.00	<input type="checkbox"/>
CCH	Crutch and Cane Holder *A1	\$60.00	<input type="checkbox"/>	CM16	Heavy Duty Cross Brace	\$95.00	<input type="checkbox"/>
ASR-A	Amputee Support - Right *A2	\$95.00	<input type="checkbox"/>	CM17	Rigid Frame	\$400.00	<input type="checkbox"/>
ASL-A	Amputee Support - Left *A2	\$95.00	<input type="checkbox"/>				
CM6-A	1-Piece Footboard	\$95.00	<input type="checkbox"/>				
SEAT, BACK & UPHOLSTERY				ADDITIONAL NOTES			
SSI	Solid Seat Pocketed Insert *A2	\$80.00	<input type="checkbox"/>				
RBCR	Raised Rigidizer	\$120.00	<input type="checkbox"/>				
TBU-SW	Adjustable Tension Back - 14" - 18"	N/C	<input type="checkbox"/>				
TBUC-SW	Adjustable Tension Back over 18" (*5)	N/C	<input type="checkbox"/>				
FOOTNOTES							
*1 - Require 20"/22" wheels		*4 - ALF/AAF are recommended		*7 - 19" seat depth by upholstery on 20" frame		*A1 - Client's Responsibility	
*2 - Require "C" arms		*5 - BCR/RBCR are recommended		*8 - Unavailable with HD2 or HD3		*A2 - Clinical Justification Req.	
*3 - 5"/6" casters recommended		*6 - CM1 Requires CM2					

Prices subject to change, an order confirmation will be sent with up-to-date pricing.

www.mapleleafwheelchair.ca