



Invacare® MyOn® Jr.
Folding Wheelchair
Suggested Canadian Price List
Effective April 7th, 2020



Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001

NAME _____ PHONE NO _____

ACCOUNT# _____ P.O. _____ DATE _____

ADDRESS _____

FRAME SELECTION

FRAME TYPE

- MYONJR MyOn Jr. Folding Wheelchair \$2,425

SEAT WIDTH RANGE

- SWR811 Adj.Seat Width 8-11" No Charge
Available with SDR810, SDR1012 and SDR1214
- SWR1215 Adj.Seat Width 12-15" No Charge
Available with SDR1113 and SDR1315

SEAT WIDTH SETTINGS

Must refer to Width, Depth and Depth Range Matrix for available options.

- W08 8" Seat Width Setting (Available w/ SWR811) No Charge
- W09 9" Seat Width Setting (Available w/ SWR811) No Charge
- W10 10" Seat Width Setting (Available w/ SWR811) No Charge
- W11 11" Seat Width Setting (Available w/ SWR811) No Charge
- W12 12" Seat Width Setting (Available w/ SWR1215) No Charge
- W13 13" Seat Width Setting (Available w/ SWR1215) No Charge
- W14 14" Seat Width Setting (Available w/ SWR1215) No Charge
- W15 15" Seat Width Setting (Available w/ SWR1215) No Charge

SEAT DEPTH RANGE

Must refer to Width, Depth and Depth Range Matrix for available options.

- SDR810 Adj.Seat Depth 8"-10" Seat Depth No Charge
Available with STABLE or DUALAMP axle position. Not available with 22" wheels, 22SP, or 22ITESPK.
- SDR1012 Adj.Seat Depth 10"-12" Seat Depth No Charge
- SDR1214 Adj.Seat Depth 12"-14" Seat Depth No Charge
- SDR1113 Adj.Seat Depth 11"-13" Seat Depth No Charge
- SDR1315 Adj.Seat Depth 13"-15" Seat Depth No Charge

SEAT DEPTH SETTINGS

Must refer to Width, Depth and Depth Range Matrix for available options.

- D08 8" Seat Depth Setting No Charge
- D09 9" Seat Depth Setting No Charge
- D10 10" Seat Depth Setting No Charge
- D11 11" Seat Depth Setting No Charge
- D12 12" Seat Depth Setting No Charge
- D13 13" Seat Depth Setting No Charge
- D14 14" Seat Depth Setting No Charge
- D15 15" Seat Depth Setting No Charge

FRONT SEAT-TO-FLOOR HEIGHT

- FSTF15 Seat-to-floor 15" No Charge
Available with 80STRSHORT0, 80STRSHORT30, 80STRLONG0, 80STRLONG30.
- FSTF16 Seat-to-floor 16" No Charge
Available with 80STRSHORT0, 80STRSHORT30, 80STRLONG0, 80STRLONG30.
- FSTF17 Seat-to-floor 17" No Charge
Not available with 18COM wheels.
- FSTF18 Seat-to-floor 18" No Charge
Not available with 18COM wheels.

REAR AXLE POSITION

- ACTIVE Active (In front of Back Cane) **Standard**
- Not Available with AA748 back type.*
- STABLE Stable (Behind Back Cane) No Charge
- DUALAMP Dual Amputee Rear Wheels Position \$114

ARM & BACK SELECTION

ARMS (PAIRS) - See next page for pictures

- TAHD T Arm Adjustable Height Desk Length w/ Locking Pins \$372
- TAHF T Arm Adjustable Height Full Length w/ Locking Pins \$372
Available with D13-D15.
- TAHJ T Arm Adjustable Height Junior \$372
- REMG Removable Mudguard \$203
- RESK Removable Skirtguard \$203
Not Available with FSTF15.
- LESSARMS Less Arms No Charge



Weight Capacity 175 lbs

MYONJR SEATING MATRIX (Canada)								
Seat Width \ Seat Depth	Seat Width							
	W08	W09	W10	W11	W12	W13	W14	W15
D08	Y	Y	Y	Y	X	X	X	X
D09	Y	Y	Y	Y	X	X	X	X
D10	Y	Y	Y	Y	X	X	X	X
D11	Y	Y	Y	Y	Y	Y	Y	Y
D12	Y	Y	Y	Y	Y	Y	Y	Y
D13	Y	Y	Y	Y	Y	Y	Y	Y
D14	Y	Y	Y	Y	Y	Y	Y	Y
D15	X	X	X	X	Y	Y	Y	Y

Y Compatible
X Not compatible


MyOn Jr Width, Depth and Depth Range Matrix									
		W8	W9	W10	W11	W12	W13	W14	W15
SDR810	D8	X	X	X	X				
	D9	X	X	X	X				
	D10	X	X	X	X				
SDR1012	D11	X	X	X	X				
	D12	X	X	X	X				
SDR1214	D12	X	X	X	X				
	D13	X	X	X	X				
SDR1113	D11					X	X	X	X
	D12					X	X	X	X
SDR1315	D13					X	X	X	X
	D14					X	X	X	X
	D15					X	X	X	X

X = Width, Depth and Depth Ranges that are available together

FSTF Height Chart				
MyOn Jr	20	22		
FSTF	5	6	5	6
15	X		X	
16	X	X	X	X
17	X	X	X	X
18	X	X	X	X
	Flush Seat Slope Only			
	Flush & 1" Seat Slope			
	Flush, 1", and 2" Seat Slope			
	1" & 2" Seat Slope			

BACK TYPE (Modular backs should only be used with straight canes)

<input type="checkbox"/> FIX057	Fixed 0° Backrest, Adj. Hgt Canes,Uph. Height 15"-17"	No Charge
<i>Modular backs should only be used with straight canes. Add 2.8" for overall cane height.</i>		
<input type="checkbox"/> FIX03P	Fixed 0° Backrest, 3 Position Uphl. Height 15"/16"/17", Cane Hgt 20"	No Charge
<i>Factory setting is at 15". Modular backs should only be used with straight canes</i>		
<input type="checkbox"/> FIX746	Fixed 7° Backrest, Adj. Height Canes, Uph. Height 14-16"	No Charge
<i>Add 2.3" for overall cane height.</i>		
<input type="checkbox"/> FIX73P	Fixed 7° Backrest, 3 Position Uphl. Height 15"/16"/17", Cane Hgt 20"	No Charge
<i>Add 2.3" for overall cane height.</i>		
<input type="checkbox"/> AA748	Adj. Angle 7° Backrest, Adj.Hgt Canes,Uph. Height 14"-18" (-6° to 16°,infinite range)	No Charge
<i>Add 2.3" for overall cane height.</i>		
<input type="checkbox"/> FOLD3P	Folding Backrest Uphl. Height 14"-16", Cane Height 20"	\$172
<input type="checkbox"/> FOLD16	Folding Bkrst Uphl. Height 14"-16" w/Adj. Hgt Push Handles	\$172
<i>Maxium Cane height will be 20" when using the 16" Upholstery height option.</i>		



CASTER & WHEEL OPTIONS

CASTER OPTIONS

<input type="checkbox"/> 5ALUM	5" Aluminum Urethane Casters	\$153
<input type="checkbox"/> AT5SOFT	5" x 1 1/2" Soft Roll Casters	No Charge
<input type="checkbox"/> AT901	5" Composite Urethane Casters	No Charge
<input type="checkbox"/> AT902	6" Composite Urethane 3 Spoke Casters	No Charge
<input type="checkbox"/> 6PNEU	6" X 1 1/4" Composite Pneumatic Casters	\$73
<input type="checkbox"/> STBKIT	Stabilizer Kit	\$83

WHEEL OPTIONS

<input type="checkbox"/> 18COM	18" Composite Tires	No Charge
<i>Select PNEU or PNEUFF in TIRE OPTIONS.</i>		
<input type="checkbox"/> 1763BK	20" Composite Pneumatic Tires	\$231
<input type="checkbox"/> 4763BK	20" Composite Pneumatic Flat Free Inserts	\$349
<input type="checkbox"/> U1763BK	20" Composite Urethane Tires	No Charge
<input type="checkbox"/> 20SP	20" Spoke	No Charge
<i>Comes standard with Pneumatic tires.</i>		
<input type="checkbox"/> 22SP	22" Spoke	No Charge
<i>Comes standard with Pneumatic tires.</i>		
<input type="checkbox"/> 1764BK	22" Composite Pneumatic Tires	\$231
<input type="checkbox"/> 4764BK	22" Composite Pneumatic Flat Free Inserts	\$349
<input type="checkbox"/> U1764BK	22" Composite Urethane Tires	No Charge
<input type="checkbox"/> 22SPITESPK	22" Spoke	No Charge
<i>Comes standard with Pneumatic tires.</i>		
<input type="checkbox"/> 22ITESPK	22" Gallop Light Spoke	No Charge
<i>Comes standard with Pneumatic tires.</i>		

TIRE OPTIONS

<input type="checkbox"/> PNEU	Pneumatic	No Charge
<i>Available with wheels: 18COM.</i>		
<input type="checkbox"/> PNEUFF	Pneumatic-Flat Free	\$128
<i>Available with wheels: 18COM.</i>		

HANDRIM OPTIONS

<input type="checkbox"/> ALUMPC	Aluminum Anodized Handrims	Standard
<i>Not available with U2222C</i>		
<input type="checkbox"/> A1512	Plastic Coated Handrims	\$114
<input type="checkbox"/> AT1225-8	Plastic Coated w/8 Vertical Projections	\$307
<i>Available on 4763BK, 1764BK, 4764BK or U1764BK</i>		
<input type="checkbox"/> L8702	Omit Handrims Only	Less net (\$15)

WHEEL LOCKS

<input type="checkbox"/> WLPULL	Pull-To-Lock Wheel Lock (Clamp-On Style)	No Charge
<i>Not available together with: FSTF15 and 22" rear wheels or FSTF16 and 24" rear wheels.</i>		
<input type="checkbox"/> WLPUSH	Push-To-Lock Wheel Lock (Clamp-On Style)	No Charge
<i>Not available with DUALAMP.</i>		
<input type="checkbox"/> WLPULLX	Pull-To-Lock Wheel Lock w/Wheel Lock Extensions (Clamp-On Style)	No Charge
<i>Not available with DUALAMP. FSTF15 and 22" rear wheels or FSTF16 and 24" rear wheels.</i>		
<input type="checkbox"/> WLPUSHX	Push-To-Lock Wheel Lock w/Wheel Lock Extensions (Clamp-On Style)	No Charge
<i>Not available with DUALAMP.</i>		



REAR AXLES

<input type="checkbox"/> PAXLE	Permanent Axle	No Charge
<i>Not available with DUALAMP.</i>		
<input type="checkbox"/> A1577	Quick Release Axles	No Charge



Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001

FRAME & UPHOLSTERY COLOUR

CONTEMPORARY

- 24P Wet Black
- 60P Silver Metallic
- 71P Silver Vein
- 105P Midnight Blue
- 115P Black Prism

TOO HOT

- 30P Sunny Yellow
- 61P Electric Red
- 119P Bubblegum Pink
- 120P Grape Madness

WAY COOL

- 62P Electric Blue
- 104P Electric Teal
- 121P Lolly Pop Blue
- 122P Cosmic Blue
- 125P Emerald Green
- 127P Grasshopper Green
- NN380 Touch-Up Paint Colours

SEAT UPHOLSTERY TYPE

- ADJWDS Adjustable Width & Depth Seat Upholstery
- L8651 Omit Seat Upholstery

BACK UPHOLSTERY TYPE

- ADJWHB Adjustable Width & Height Nylon Type
Includes backbrace. Not available with stroller handle.
- L8601 Omit Back Upholstery
Does not include backbrace.

POSITIONING STRAPS

SEAT POSITIONING *(Available on W11 and wider.)*

- A1515 Auto Style Seat Positioning Strap \$84
- A1321 Hook & Loop Seat Positioning Strap \$66

CHEST POSITIONING *(Available on W11 - W15.)*

- 9515 Auto Style Chest Positioning Strap \$84
- 7321 Hook & Loop Chest Positioning Strap \$66

MISCELLANEOUS ACCESSORIES

- A1580P Anti-Tippers w/ Tipping Aid (pair) \$114
- P1580 Rear Anti Tippers (pair) \$109
- A1569 Crutch and Cane Holder \$66
- ACTRY Acrylic Tray *(fits 12"-14" Widths)* \$270

Not available with LESSARMS.

- ACTRSA Right Side Swing Away Acrylic Tray *(fits 8"-15" widths)* \$135
- ACTLSA Left Side Swing Away Acrylic Tray *(fits 8"-15" widths)* \$135
- A1307SM Spoke Protectors - Smiley Face Design \$207

Available with 20SP, 22SP, 22ITESPK.

- A1307SP Spoke Protectors - Spider Web Design \$207

Available with 20SP, 22SP, 22ITESPK.

- A1307FL Spoke Protectors - Floral Design \$207

Available with 20SP, 22SP, 22ITESPK.

- STRHNDL Stroller Handle *(clamp on)* \$156

Not available with upholstery back.

- MYCVR Monkey Cover for 90° Straight Lift Off Hanger \$19

Available on 90STRLIFT

- BRCVR Bear Cover for 90° Straight Lift Off Hanger \$19

Available on 90STRLIFT

FRONT RIGGING OPTIONS

FRONT RIGGING (PAIRS)

- 90STRLIFT 90° Straight Lift Off Hanger \$150
- 80STRSHORT0 80° Straight Short Swing In/Out Hng, 0° Std Position *(blk frm colour)* \$150
- 80STRSHORT30 80° Straight Short Swing In/Out Hng, 30° Abduction Pos. *(blk frm colour)* \$150
- 80STRLONG0 80° Straight Long Swing In/Out Hng, 0° Std Position \$150
- 80STRLONG30 80° Straight Long Swing In/Out Hng, 30° Abduction Position \$150
- ELVOFF Elevating Legrest Lift Off *(0° - 80°, includes calf pads)* \$416

Available on W13-W15.

- ELVOFFWO Elevating Legrest Lift Off *(0° - 80° w/o calf pads, includes calf strap)* \$416

Available on W12-W15.

- NO FRONTRIG No Front Rigging Selection No Charge

FOOTPLATE TYPE (PAIRS)

- FLIPFTS Flip Up Footplates No Charge

Available on 12"-15" widths. Not Available on 90STRLIFT.

- FLPDADJ Flip Up Depth & Angle Adj. Footplates \$135

Available on 8" or 9" widths.

- 1PCFTPT One Piece Footplate \$135

Not available on ELVOFF or ELVOFFWO.

Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001





No Charge
No Charge
No Charge
No Charge
No Charge

No Charge
No Charge
No Charge
No Charge

No Charge
No Charge
No Charge
No Charge
No Charge
\$44

Standard
Less net (\$15)

Standard
Less net (\$20)



FOOTREST/LEGREST ACCESSORIES

- A1337 Calf/leg Strap \$63
- A9342 "H" Heel Strap \$94
- Available on W12-W15.*
- HLOOP Heel Loops \$79
- Available on W11 or greater.*
- FRMPRBK ChairArmor Black Neoprene Frame Protectors No Charge
- Fits on 80° and 90° Front rigging. Not Available on Elevating Legrests.*


MyOn Jr STF Chart - general guideline

Flush Seat Slope		Wheel size and Axle position		Fork Position & Caster Size F1=1 hole fork, F2=2 holes, F3=3 holes			
Front STF	Rear STF	20"	22"	Fork in "UP" position		Fork in "DOWN" position	
15" FSTF	15" RSTF	axle 2	axle 1	F1			
16" FSTF	16" RSTF	axle 3	axle 2			F1	
17" FSTF	17" RSTF	axle 4	axle 3	F2 bottom	F2 top		
18" FSTF	18" RSTF	axle 5	axle 4		F2 bottom	F2 bottom	F2 top

1" Seat Slope		Wheel size and Axle position		Fork Position and Caster Size			
Front STF	Rear STF	20"	22"	Fork in "UP" position		Fork in "DOWN" position	
15" FSTF	14" RSTF	axle 1		F1			
16" FSTF	15" RSTF	axle 2	axle 1			F1	
17" FSTF	16" RSTF	axle 3	axle 2	F2 bottom	F2 top		
18" FSTF	17" RSTF	axle 4	axle 3		F2 bottom	F2 bottom	F2 top

2" Seat Slope		Wheel size and Axle position		Fork Position and Caster Size			
Front STF	Rear STF	20"	22"	Fork in "UP" position		Fork in "DOWN" position	
15" FSTF	13" RSTF						
16" FSTF	14" RSTF	axle 1				F1	
17" FSTF	15" RSTF	axle 2	axle 1	F2 bottom	F2 top		
18" FSTF	16" RSTF	axle 3	axle 2		F2 bottom	F2 bottom	F2 top

			Pivot Tube part number	Knee-to-heel range (from seat rail to top of foot plate)
MYONJR	80STRSHRT0 , 80STRSHRT30, 80STRLONG0, 80STRLONG30	FLIPFTS (Composites), FLPDADJ (Angle/Depth Adjustable), and 1PCFTPT (one piece footplate)	Full range of seat heights with tube 1511849	80 deg. Short hanger with FLIPFTS= 7.5 - 12.75"
				80 deg. Short hanger with FLPDADJ= 7.5 - 12.75"
				80 deg. Long hanger with FLIPFTS= 12" - 15"
				80 deg. Long hanger with FLPADJ= 12" - 15"
90STRLIFT	FLPDADJ 1PCFTPT		1506244-0699 for all lengths	5.25" - 10.75"
				Same tube for height range; tube depends on width (comes with plate)
				5.25" - 10.75"



Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
WWW.INVACARE.CA