



REAR FLOOR TO SEAT HEIGHT TABLE

REAR WHEEL DIAMETER XL5 <i>select desired height</i>		
POSITION	20"	22"
A	13 <input type="checkbox"/>	14 <input type="checkbox"/>
B	13 ½ <input type="checkbox"/>	14 ½ <input type="checkbox"/>
C	14 <input type="checkbox"/>	15 <input type="checkbox"/>
D	14 ½ <input type="checkbox"/>	15 ½ <input type="checkbox"/>
E	15 <input type="checkbox"/>	16 <input type="checkbox"/>
F	15 ½ <input type="checkbox"/>	16 ½ <input type="checkbox"/>
G	16 <input type="checkbox"/>	17 <input type="checkbox"/>
H	16 ½ <input type="checkbox"/>	17 ½ <input type="checkbox"/>
I	17 <input type="checkbox"/>	18 <input type="checkbox"/>
J	17 ½ <input type="checkbox"/>	18 ½ <input type="checkbox"/>



Note - The measures are for reference purposes only. They depend on the wheel diameter, tilt angle and the rear wheel axle location.

REFERENCE SCHEME





A: SEAT-TO-FLOOR HEIGHT : _____

B: REAR SEAT-TO-FLOOR HEIGHT: _____

REAR WHEEL POSITION

Choose the rear wheel desired position



NOTES

AADL EXCLUSION LIST

PEDIATRIC XL5 - NOT FUNDED OPTION LIST

- Fixed Extension Handle
- Plastic Coated Handrim
- Center-Pull Auto Style Hip Belt