



Squiggles Stander - AADL

October 2020



Mark For:	Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order
Date:	PO#:
Dealer Acct #:	ADDITIONAL SHIPPING INFORMATION
Dealer:	Ship To:
Dealer Contact:	Attention:
Dealer Address:	Address:
Dealer City: PV: PC:	Address:
Dealer Phone: Fax:	Ship To City: PV: PC:
Confirmation Email:	Ship To Phone: Fax:
Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	

Squiggles Stander

Weight Capacity: 48.4 lbs

SECTION A - PACKAGES

Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps, and Product Assembly. Substitutions are not available for Packages.

	Item Code	Description	Retail
<input type="checkbox"/>	PSSTB	Squiggles Stander Package - Blue	\$5,060.00
<input type="checkbox"/>	PSSTBV	Squiggles Stander Package - Blue Vinyl	\$5,060.00
<input type="checkbox"/>	PSSTG	Squiggles Stander Package - Green	\$5,060.00
<input type="checkbox"/>	PSSTO	Squiggles Stander Package - Orange	\$5,060.00
<input type="checkbox"/>	PSSTP	Squiggles Stander Package - Pink	\$5,060.00

SECTION B - MADE-TO-ORDER (Required if you did not choose a package in Section A)

Standard Model

Selection Required. Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.

	Item Code	Description	Retail
<input type="checkbox"/>	118-600	Standard Model with Support Shell	\$1,995.00
	118-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning	\$1,484.00

Positioning Pads

Selection Required. Includes chest and hip pad covers and knee support covers.

	Item Code	Description	Retail
<input type="checkbox"/>	118-762	Squiggles Stander Positioning Pads - Blue	\$562.00
<input type="checkbox"/>	118-762-03VIN	Squiggles Stander Positioning Pads - Blue Vinyl	\$562.00
<input type="checkbox"/>	118-760	Squiggles Stander Positioning Pads - Green	\$562.00
<input type="checkbox"/>	118-761	Squiggles Stander Positioning Pads - Orange	\$562.00
<input type="checkbox"/>	118-763	Squiggles Stander Positioning Pads - Pink	\$562.00



Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001

Customer Service: (800) 263-3390 Fax: (800) 561-5834 E-mail: cscanada@sunmed.com www.sunrisemedical.ca

© 2020 Sunrise Medical Canada Inc. MK-129858-CAN Rev (1).

Accessories					
	Item Code	Description	Quantity:		Retail
<input type="checkbox"/>	102-1640	Size 1 Knee Cup		(1 or 2)	\$111.00
<input type="checkbox"/>	102-3640	Size 3 Knee Cup		(1 or 2)	\$111.00
<input type="checkbox"/>	118-619	Tray - Gray			\$374.00
<input type="checkbox"/>	120-270	Padded Tray Insert - Gray			\$305.00
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray			\$281.00
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)			\$210.00
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)			\$210.00
<input type="checkbox"/>	120-816	Footplate Raiser 1" (Each)		(1 to 3)	\$84.00

Headrest

Cushion selection required with hardware

	Item Code	Description	Retail
<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware	\$340.00
<input type="checkbox"/>	137-800-03	Contoured Headrest Cushion - Blue	\$167.00
<input type="checkbox"/>	137-800-03VIN	Contoured Headrest Cushion - Blue Vinyl	\$167.00
<input type="checkbox"/>	137-800-01	Contoured Headrest Cushion - Green	\$167.00
<input type="checkbox"/>	137-800-02	Contoured Headrest Cushion - Orange	\$167.00
<input type="checkbox"/>	137-800-04	Contoured Headrest Cushion - Pink	\$167.00
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	\$228.00
<input type="checkbox"/>	117-850-03	Flat Headrest Cushion - Blue	\$114.00
<input type="checkbox"/>	117-850-03VIN	Flat Headrest Cushion - Blue Vinyl	\$114.00
<input type="checkbox"/>	117-850-01	Flat Headrest Cushion - Green	\$114.00
<input type="checkbox"/>	117-850-02	Flat Headrest Cushion - Orange	\$114.00
<input type="checkbox"/>	117-850-04	Flat Headrest Cushion - Pink	\$114.00
<input type="checkbox"/>	118-667	Lateral Bracket Assembly Hardware (Pair) <i>Requires Flat Headrest Support Hardware and Cushion</i>	\$162.00
<input type="checkbox"/>	137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	\$298.00
<input type="checkbox"/>	137-807-03VIN	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)	\$298.00
<input type="checkbox"/>	137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	\$298.00
<input type="checkbox"/>	137-807-02	Flat Headrest Lateral Cushions - Orange (Pair)	\$298.00
<input type="checkbox"/>	137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	\$298.00



Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001

Customer Service: (800) 263-3390 Fax: (800) 561-5834 E-mail: cscanada@sunmed.com www.sunrisemedical.ca

© 2020 Sunrise Medical Canada Inc. MK-129858-CAN Rev (1).

Technical Specifications

Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kgs
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5" Size 3: 3.5"	Size 1: 63 mm Size 3: 89 mm



Select store location - AADL Numbers

- Calgary - AADL 00000 78000 003
- Edmonton - AADL 00000 78000 001
- Grande Prairie - AADL 0000 665538 001



Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3 • Concord, ON • Canada L4K 4V3