

**VGP VOYAGER PLUS Wheelchair (Manual Wheelchair)**



Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_

FRAME COLOUR		Select one below	1
<input type="checkbox"/>	VOYAGER PLUS FRAME - Black		1450

SEAT WIDTH		Select one below	2
<input type="checkbox"/>	SW14"-20"	Specify _____	N/C

SW14"-20" Light Weight Aluminum Frame

<input type="checkbox"/>	SW21"-28"	Specify _____	N/C
--------------------------	-----------	---------------	-----

SW21"-28" Must Include Plus Package 1, 2 Or 3

SEAT DEPTH		Select one below	3
------------	--	------------------	---

<input type="checkbox"/>	SDF18"	(By Frame)	N/C
<input type="checkbox"/>	SD16"-18"	Specify _____ (By Upholstery)	N/C
<input type="checkbox"/>	SDF20"	(By Frame)	350
<input type="checkbox"/>	SD18"-20"	Specify _____ (By Upholstery)	N/C
<input type="checkbox"/>	SDF22"	(By Frame)	350
<input type="checkbox"/>	SD20"-22"	Specify _____ (By Upholstery)	N/C

BACK		Select one below	4
------	--	------------------	---

<input type="checkbox"/>	BH15"-20"	Specify _____	N/C
<input type="checkbox"/>	BH21"+	Specify _____	120
<input type="checkbox"/>	AB	Angle Adj. Back	N/C
<input type="checkbox"/>	ENDO*	Endomorph Canes W ( )	N/C

Specify Width

Endomorph Canes require TA or TAS Arm Rest

STROLLER HANDLE		Select one below	5
-----------------	--	------------------	---

<input type="checkbox"/>	SHS	Removable Stroller Handle	N/C
<input type="checkbox"/>	SHR	Removable Raised Stroller Handle	N/C

FOOT REST		Select one below	6
-----------	--	------------------	---

<input type="checkbox"/>	F6	60 Degree Hanger	N/C
<input type="checkbox"/>	F7	70 Degree Hanger	N/C
<input type="checkbox"/>	F8	80 Degree Hanger	N/C
<input type="checkbox"/>	EA	Elevating Legrest (A2)	404
<input type="checkbox"/>	AL	Articulating Legrest (A2)	450

FOOT REST AND FOOT PLATES		Select one below	7
---------------------------	--	------------------	---

<input type="checkbox"/>	PPL	Plastic	N/C
<input type="checkbox"/>	PAF	Angle Adjustable - Extra Large	150
<input type="checkbox"/>	PXL	Extra Large - Alum (F6A,F7A,F8A)	150
<input type="checkbox"/>	1FB	One Piece Foot Board	80

ARM REST		Select one below	8
----------	--	------------------	---

<input type="checkbox"/>	FA	Flip Back Ht. Adj./ Removable	240
<input type="checkbox"/>	TA*	T-Arm Ht. Adj./ Removable	240
<input type="checkbox"/>	FAS	Flip Back (Short) Ht. Adj./ Remove	240
<input type="checkbox"/>	TAS*	T-Arm (Short) Ht. Adj./ Removable	240
<input type="checkbox"/>	GWAPS	Gel Waterfall Arm Pads (A1)	150

14" - 20" Wide Light Weight Aluminum Frame Only

Please Note : A1 = Clients Responsibility

Please Note : A2 = Clinical Justification

**SPECIFICATIONS**

- LIGHT WEIGHT ALUMINUM FRAME AND CROSS BRACES
- STANDARD 18" BACK HEIGHT
- 20", 22", OR 24" NYLON MOULDED REAR WHEELS
- EXTRUDED ALUMINUM FORKS
- CHROME HAND RIMS
- 5", 6", OR 8" CASTERS
- BLACK NYLON UPHOLSTERY
- POWDER COATED FRAME
- FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
- DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
- PRECISION BEARINGS
- MOULDED PLASTIC FOOTPLATES
- SEAT WIDTH 14" TO 18"
- BLACK PADDED ARM PADS
- DOUBLE CROSS BRACE

**SEAT TO FLOOR SPECIFICATIONS**

**(LOW SEAT SETTING)**

	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	SL 12.25	20"	5"	XL SHORT
<input type="checkbox"/>	SL 13.00"	20"	5"	XL SHORT
<input type="checkbox"/>	LL 13.75"	20"	5"	XL SHORT
<input type="checkbox"/>	L 14.50"	20"	6"	SHORT
<input type="checkbox"/>	L 14.50"	22"	5"	SHORT
<input type="checkbox"/>	L14.50"	22"	6"	SHORT
<input type="checkbox"/>	LH 15.25"	20"	5"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	20"	6"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	5"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	6"	SHT/LONG
<input type="checkbox"/>	SH 16.00"	22"	6"	LONG

**(STANDARD SEAT SETTING)**

<input type="checkbox"/>	MH 16.75	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	22"	8"	LONG
<input type="checkbox"/>	H 17.50"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	24"	8"	LONG
<input type="checkbox"/>	TH18.25"	24"	8"	LONG
<input type="checkbox"/>	TT 19.00"*	24"	8"	LONG

\*TT 19.00' requires taller caster journal

REAR WHEELS			Select one below	9
<input type="checkbox"/>	RWT0	20" MAG Solid Treaded Tire	N/C	
<input type="checkbox"/>	RWT2	22" MAG Solid Treaded Tire	N/C	
<input type="checkbox"/>	RWT4	24" MAG Solid Treaded Tire	N/C	
<input type="checkbox"/>	RWP2	22" Pyramid Tire	80	
<input type="checkbox"/>	RWP4	24" Pyramid Tire	80	
<input type="checkbox"/>	RWA2	22" Air Tire	120	
<input type="checkbox"/>	RWA4	24" Air Tire	120	
<input type="checkbox"/>	AQR	Axles – Quick Release (349 Lbs)	N/C	
<input type="checkbox"/>	PCHR	Plastic Coated Hand Rims	N/C	

14" - 20" Wide Light Weight Aluminium Frame Only

CASTER			Select one below	10
<input type="checkbox"/>	FC1A	5" X 1.25" Caster	N/C	
<input type="checkbox"/>	FC1B	6" X 1.25" Caster	N/C	
<input type="checkbox"/>	FC1C	8" X 1.25" Caster	N/C	
<input type="checkbox"/>	FC2A	5" X 1.75" Caster	80	
<input type="checkbox"/>	FC2B	6" X 1.75" Caster	80	
<input type="checkbox"/>	FC2C	8" X 1.75" Caster	80	
<input type="checkbox"/>	FC3A	5" X 1.25" Air Caster	80	
<input type="checkbox"/>	FC3B	6" X 1.25" Air Caster	80	
<input type="checkbox"/>	FC3C	8" X 1.25" Air Caster	80	

WHEEL LOCK			Select one below	11
<input type="checkbox"/>	ALPL	Aluminium Pull To Lock	N/C	
<input type="checkbox"/>	ALPH	Aluminium Push To Lock	N/C	
<input type="checkbox"/>	STPL	Steel Pull To Lock	N/C	
<input type="checkbox"/>	STPH	Steel Push To Lock	N/C	
<input type="checkbox"/>	BRXA	6" Brake Extension	40	
<input type="checkbox"/>	BRXB	9" Brake Extension	60	

SEATING - BACKS			Select one below	12
			W ( ) x H ( )	
<input type="checkbox"/>	TAB	Tension Adjustable Back		180

ACCESSORIES				14
<input type="checkbox"/>	SBA	Seat Belt (Auto Style)		22
<input type="checkbox"/>	SBB	Seat Belt (Velcro Style)		15.75
<input type="checkbox"/>	CNC	Crutch And Cane Holder (A1)		60
<input type="checkbox"/>	IVA	Adjustable I V Pole (A1)		130
<input type="checkbox"/>	NLSP	Nylon Leg Strap Padded (A1)		50
<input type="checkbox"/>	VGPAMPL	Amputee Support Left		45
<input type="checkbox"/>	VGPAMPR	Amputee Support Right		45
<input type="checkbox"/>	VGPOXY	Oxygen Tank Holder		130
<input type="checkbox"/>	VGPATR	Anti Tippers		88

Please Note : A1 = Clients Responsibility

PLUS PACKAGES		Select one below	15
<input type="checkbox"/>	<b>PLUS PACKAGE 1: AADL CODE: W101</b>		
<b>CLIENT 250 TO 349 LBS - \$880.00</b>			
INCLUDES THE FOLLOWING:			
STEEL FRAME, DOUBLE CROSS BRACES WITH HEAVY DUTY BACK AND SEAT UPHOLSTERY			
SEAT WIDTH INCREASE FROM 16" TO 24"			
SEAT DEPTH INCREASE FROM 16" TO 22"			
5,6 OR 8 X 1 1/4" HD CASTERS			
*** WELDED CASTER JOURNAL ***			
<input type="checkbox"/>	<b>PLUS PACKAGE 2: AADL CODE: W102</b>		
<b>CLIENT 350 TO 449 LBS - \$1140.00</b>			
INCLUDES THE FOLLOWING:			
STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY			
SEAT WIDTH INCREASE FROM 18" TO 26"			
SEAT DEPTH INCREASE FROM 16" TO 22"			
5,6 OR 8 X 1 1/4" HD CASTERS			
HD STEEL FORKS			
*** WELDED CASTER JOURNAL ***			
<input type="checkbox"/>	<b>PLUS PACKAGE 3: AADL CODE: W103</b>		
<b>CLIENT 450 TO 600 LBS - \$1300.00</b>			
INCLUDES THE FOLLOWING:			
STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY			
SEAT WIDTH INCREASE FROM 20" TO 28"			
SEAT DEPTH INCREASE FROM 16" TO 22"			
5,6 OR 8 X 1 1/4" HD CASTERS			
HD STEEL FORKS			
*** WELDED CASTER JOURNAL ***			
<input type="checkbox"/>	FIXED REINFORCED ONE PIECE BACK POST	\$260.00	
(SPECIFY ANGLE _____)			
<input type="checkbox"/>	HEAVY DUTY TENSION ADJ. BACK UPH.	\$220.00	

Please Note : Reference Page 12 of the AADL Wheelchair APL for options and catalogue numbers for A wheelchairs.

Notes: