## ARJOHUNTLEIGH GETINGE GROUP





...with people in mind

## PORTABLE CEILING LIFT VOYAGER


The Voyager® portable ceiling lift combines a solidly engineered design with features you would never expect to find in a portable lift.



Our unique Quick Release feature allows the user to disengage the strap and guckly install or remove the lift from the track.



The handset offers a degree of protection against splashes (IPX4)



A 24-volt battery allows you to perform all your tasks worry-free. When necessary, it recharges in less than 2 hours.



Combined with the Quick Release function, the arm extension gives the user the ability to detach the lift from the track. At the same time, unnecessary strain is avoided for the caregiver.

## Compatible Slings - ArjoHuntleigh Ergofit slings are designed to be easily installed and removed.



Hammock Sling THA model

Weight

**Product Information** 

Soft-start & soft-stop motor control **Emergency Lowering Device** 

Lifting Capacity

**Emergency Stop Current Limiter** 



Quick-Fit Sling TIR model



Hygienic Sling THY model



Walking Sling	
TEM model	

1.7 in/sec (4.3 cm/sec)
88" (2.2 m)

Use only ArjoHuntleigh designed parts, i.e. parts specifically designed for the purpose, on equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development, we reserve the right to modify designs and specifications without prior notice

12 lb (5.5 kg)

440 lb (200 kg)

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies. © ArjoHuntleigh, 2014

## **ArjoHuntleigh**