

LTC 4000 Ultra 55 Mattress

Ultra Stable Firm Memory Foam

Firm viscoelastic foam plus a high resiliency transitional layer with full height transfer borders ensure stability during transfers and when performing ADL's. Helps prevent entrapment between mattress and bedside rails.

2" Medium-plush Viscoelastic Foam Comfort Layer

2" High Resiliency Transitional Layer

2" Firm Core

Patient Considerations for this Mattress:

Pressure Distribution	> E
Firmness	> F
Temperature	>(
Transfers	> E
Bed Mobility	> \
Other	>(

- > Best choice for individuals at risk for pressure sore development
- Recommended for individuals requiring a firm mattress
- Good choice for prevention of heat build-up and sweating
- > Best choice for bedside stability, preventing entrapment
- > Viscoelastic top layer may impair bed mobility for some individuals
- > Good choice for heavier individuals or individuals preferring a firm mattress

(1) (1)

Specifications:

Top Foam Layer	2" medium-plush viscoelastic foam comfort layer plus 2" high resiliency transitional layer
Core Foam	> Firm core
Transfer Borders	> Full height exterior foam adds firmness along edges
Cover	 Fluid-proof, non-allergenic, antibacterial, antimicrobial cover Also available with organic cotton cover (not fluid-proof)
Misc	> Designed for use on all hospital beds.

Standard Sizes: (Custom sizes also available by quotation.)

- 80" × 35" × 6" • 80" × 36" × 6" • 80" × 54" × 6"
- 80" x 39" x 6" 80" x 59" x 6"

Recommended Weight Capacity:

450 lbs

100% Customer Satisfaction Guarantee:

The warranty period is five years; two years full plus three years pro-rated.

This information provides general guidance. EMP recommends specialized sleep surfaces be selected with the help of a qualified medical professional.