

SPECIFICATIONS

	Go-Go® Ultra X	Go-Go® LX w/ CTS Suspension	Go-Go® Elite Traveller®	Go-Go® Sport	Go-Go® Folding Scooter
Max. weight capacity:	260 lbs.	300 lbs.	300 lbs.	325 lbs.	250 lbs.
Maximum speed: ¹	Up to 4 mph	Up to 4 mph	Up to 4 mph	Up to 4.7 mph	Up to 4 mph
Ground clearance: ²	2.25" at midframe; 4.3" at motor	1.5" at motor	1.38" at motor	3" at midframe; 1.5" at motor	2.3"
Turning radius: ²	3W: 33" / 4W: 44.125"	3W: 33.5" / 4W: 44.75"	3W: 33" / 4W: 44"	3W: 36" / 4W: 50.5"	32"
Overall length: ²	3W: 37" / 4W: 39.5"	3W: 39.75" / 4W: 41.75"	3W: 37" / 4W: 39.5"	3W: 42" / 4W: 43"	37"
Overall width: ²	19.5"	20.5"	19.5"	21.25"	19"
Battery charger:	Off-board, 2 amp	Off-board, 3.5 amp	Off-board, 3.5 amp	Off-board, 3.5 amp	2-amp, off-board (SLA) or 2.5 amp, off-board (lithium ion)
Per-charge range: ^{1,5}	3W: Up to 6.9 miles 4W: Up to 7.2 miles	3W: Up to 11.7 miles 4W: Up to 14.6 miles	3W: Up to 14.6 miles 4W: Up to 13.9 miles	3W: Up to 16 miles 4W: Up to 16.2 miles	Up to 13 miles (lithium ion) Up to 9.3 miles (SLA)
Battery weight: ⁴	9.3 lbs. each ⁶	13.3 lbs. each ⁶	13.3 lbs. each ⁶	13.3 lbs. each ⁶	12 AH SLA - 9 lbs. ea. 1.5 AH lithium ion - 5 lbs. ea.
Batteries (2 required): ³	12 AH Interceptor	22 AH Interceptor	22 AH Interceptor	22 AH Interceptor	12 AH SLA (2 required) 11.5 AH lithium ion (1 -required)
Weight w/o battery pack: ²	3-wheel: 76 lbs.; 4-wheel: 80 lbs.	3-wheel: 85 lbs.; 4-wheel: 87 lbs.	3-wheel: 76 lbs.; 4-wheel: 80 lbs.	3-wheel: 93 lbs.; 4-wheel: 99 lbs.	51 lbs.
Heaviest piece: ²	3-wheel: 29.5 lbs. (front section) 4-wheel: 33 lbs. (front section)	3-wheel: 34 lbs. (front section) 4-wheel: 37 lbs. (front section)	3-wheel: 29 lbs. (front section) 4-wheel: 34 lbs. (front section)	3-wheel: 35 lbs. (front section) 4-wheel: 44 lbs. (front section)	51 lbs. (unit does not disassemble -unit weight without battery pack)



- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.
- AGM or gel-cell type recommended.
- Battery weight may vary based on manufacturer.
- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

Dealer Information

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

CA-GOGOTRI_10-14-19

GOGO
TRAVEL MOBILITY



Get Out
There!



pridemobility.ca

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

GO GO[®] ULTRA X



Color



Red

ALSO AVAILABLE
IN 4-WHEEL

FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- A standard front basket provides ample storage
- Auto-latching lockup mechanism

GO GO[®] LX WITH CTS SUSPENSION



ALSO AVAILABLE
IN 3-WHEEL

FEATURES

- CTS Suspension (Comfort-Trac front and rear suspension)
- Feather-touch disassembly permits simple frame separation with only one hand
- Tiller-mounted LED lighting safely illuminates your path
- Delta tiller with wrap-around handles, for left or right-handed operation, lets you rest your wrist comfortably while driving

Included Colors



Red Blue

Optional Colors



Orion Yellow Onyx Black Pearl Pink Sunburst Orange

GO GO[®] ELITE TRAVELLER[®]



NEW DELTA
TILLER DESIGN!

ALSO AVAILABLE
IN 4-WHEEL

GO GO[®] SPORT



ALSO AVAILABLE
IN 3-WHEEL

FEATURES (Go-Go Elite Traveller and Go-Go Sport)

- Feather-touch disassembly permits simple frame separation with only one hand
- Tiller-mounted LED lighting safely illuminates your path
- Delta tiller with wrap-around handles, for left or right-handed operation, lets you rest your wrist comfortably while driving
- Charger port located on the tiller enables you to conveniently charge your scooter
- A standard front basket provides ample storage

Included Colors



Red Blue

GO GO[®] FOLDING SCOOTER



Color



White

FEATURES

- Unique folding design
- Folds in just three easy steps
- Height-adjustable tiller
- Easy to transport at home, hotels, or airports
- Can be stowed upright or laid down

To Fold:



Fold Seat



Pull Up Lever to Fold



Fold and Lock Tiller

To Unfold:



Pull Lever



Unfold Scooter



Unfold Seat; Lock in Place