

Aluminum Rehab Shower Commode Chair with Four Rear-locking Casters

HCPCS:



FIGURE A

- Can be positioned over a standard toilet or used as a portable self-contained commode
- Floor to seat height can be adjusted 2" without tools
- White aluminum frame is attractive and rust resistant
- Padded seat back and arms
- Height adjustable
- Comes with swing-away footrests

PRODUCTS	Description
NRS185007	Shower&Commode Chair,1/ea

SPECIFICATIONS	Description
Back Height	19 "
Carton Shipping Weight	40 lbs.
Overall Height	39" - 41"
Overall Width Open	22.25"
Seat Dimensions	18" (W) x 17" (D)
Seat To Floor Height	20" - 22"
Warranty	Limited Lifetime
Weight	35.25 lbs.
Weight Capacity	275 lbs

PARTS	Description
STDS851-SA	Velcro Seat Belt Self Attaching,1ea, EA
11106-1	Commode Bucket,7.5qt.1/ea, EA
NRS185006-21	BackLockingKnobForNRS185006/7, 1/pr, PR
NRS185006-14	DowntubePlasticInsrtsForNRS185006/7,1/ea, EA
NRS185006-12	FootplateForNRS185006/7, 1/pr, PR
NRS185006-16	FootrstPivotSaddleW/BoltNRS185006/7,1/ea, EA
NRS185006-22	AntiRattleInsertForNRS185006/7, 1/ea, EA
NRS185006-17	FootrestFrameLatchForNRS185006/7, 1/ea, EA
NRS185006-15	FootrstHgtAdjKnobW/BoltNRS185006/7,1/ea, EA
NRS185006-13	AngleAdjFtplateMntW/HdwrNRS185006/7,1/ea, EA
NRS185006-11	SwingawayFootrestForNRS185006/7, 1/pr, PR
NRS185006-09	FrameClipForNRS185006/7, 4/pkg, PKG
NRS185006-08	CasterWithLegForNRS185006/7, 1/ea, EA
NRS185006-05	ArmrestBracketLeftForNRS185006/7, 1/ea, EA

PARTS**Description**

NRS185006-06	ArmrestBracketRightForNRS185006/7, 1/ea, EA
NRS185006-03	BackInsertForNRS185006/7, 1/ea, EA
NRS185006-04	ArmpadForNRS185006/7, 1/ea, PR
NRS185006-02	SeatClipW/ScrewsForNRS185006/7, 1/ea, EA
NRS185006-01	SeatWithInsertForNRS185006/7, 1/ea, EA