

SEAT WIDTH RANGE SETTINGS (continued)

Note: Frame widths measured from outside of seat rails.

- NXT-30260 16" / 17" / 18" Width Range No Charge
- NXT-W16 Set at 16" Wide
- NXT-W17 Set at 17" Wide
- NXT-W18 Set at 18" Wide
- NXT-30261 17" / 18" Width Range No Charge
- NXT-W17 Set at 17" Wide
- NXT-W18 Set at 18" Wide
- NXT-30262 18" - Not width adjustable No Charge
- NXT-W18 Set at 18" Wide

SEAT DEPTH RANGE SETTINGS

Note: NXT has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating / hardware thickness (typically 2") to client's actual seat depth.

Example: seat depth of 12" + back / hardware thickness of 2" = back post location set at 14"

Note: Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.

- NXT-30511 **SHORT FRAME & SHORT TILT PIVOT BRACKET** No Charge
- NXT-30504 Set at 12" Depth No Charge
- NXT-30505 Set at 13" Depth No Charge
- NXT-30506 Set at 14" Depth No Charge
- NXT-30521 Set at 15" Depth No Charge
- NXT-30522 Set at 16" Depth No Charge
- NXT-30523 17" - Must select 2" extension No Charge
- NXT-30524 18" - Must select 2" extension No Charge
- NXT-30525 2" Extension - Pair \$33

Note: Floor to seat frame height range is 12.5 to 16.5".

Note: Seat frame widths of 19" and 20" not available with short frame.

Note: 22" & 24" wheels not available with short frame and short tilt bracket.

- NXT-30512 **MEDIUM FRAME & SHORT TILT PIVOT BRACKET** No Charge
- NXT-30506 14" No Charge
- NXT-30521 15" No Charge
- NXT-30522 16" No Charge
- NXT-30523 17" No Charge
- NXT-30524 18" No Charge
- NXT-30527 19" - Must select 2" extension No Charge
- NXT-30528 20" - Must select 2" extension \$625
- NXT-30529 21" - Must select 3" extension \$625
- NXT-30525 2" Extension - Pair \$33
- NXT-30530 3" Extension - Pair \$33

Note: Floor to seat frame height range is 12.5 to 16.5".

Note: 22" & 24" wheels not available with medium frame and short tilt bracket.

- NXT-30540 **MEDIUM FRAME & TALL TILT PIVOT BRACKET** No Charge
- NXT-30506TP 14" No Charge
- NXT-30521TP 15" No Charge
- NXT-30522TP 16" No Charge
- NXT-30523TP 17" No Charge
- NXT-30524TP 18" No Charge
- NXT-30527TP 19" - Must select 2" extension No Charge
- NXT-30528TP 20" - Must select 2" extension \$625
- NXT-30529TP 21" - Must select 3" extension \$625
- NXT-30535TP 22" - Must select 4" extension \$625
- NXT-30525 2" Extension - Pair \$33
- NXT-30530 3" Extension - Pair \$33
- NXT-30531 4" Extension - Pair \$33

Note: Floor to seat frame height range is 14 to 18".

Note: Limits posterior tilt to 40°, includes split collars for adjustment.

Note: Medium with Tall Tilt and Long frames only.

Note: NXT has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame.

For back post location add seating / hardware thickness (typically 2") to client's actual seat depth.

Example: seat depth of 12" + back / hardware thickness of 2" = back post location set at 14"

Seat to Floor Height Requested _____ " *REQUIRED*

PUSH HANDLES

- NXT-30518 Standard 90° Bend Push Handles and Back Post Tubes No Charge
- NXT-30519 Ergo Bend Push Handles and Back Post Tubes No Charge
- NXT-30532 Seating Dynamics Rocker Back \$2,168

Note: Angle adjustable from 72° to 108° in 6° increments.

Note: Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.

Note: Locks into a non-dynamic position when necessary.

Note: Does not allow the back post tubes to fold forward.

BACK POST HEIGHT

Note: Back Posts are measured from top of seat rail to top of handle bend.

- NXT-30515 Short Push Handle - 20" high No Charge
- NXT-30514 Medium Push Handle - 23" high No Charge
- NXT-30516 Tall Push Handle - 26" high No Charge

NXT Wheel / Caster Seat to Floor Heights							
Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5"	18"	5"	Tall	16.5"
12"	3"	Tall	14"	18"	5"	Short	15.5"
12"	3"	Short	13"	18"	5"	Tall	17"
12"	3"	Tall	14.5"	18"	6"	Short	15.5"
12"	4"	Short	13.5"	18"	6"	Tall	17"
12"	4"	Tall	15"	20"	4"	Short	14"
12"	5"	Short	15.5"	20"	4"	Tall	15.5"
12"	5"	Tall	17"	20"	4"	Short	14.5"
12"	6"	Short	16"	20"	4"	Tall	16"
12"	6"	Tall	17.5"	20"	4"	Short	15"
16"	3"	Short	12.5"	20"	4"	Tall	16.5"
16"	3"	Tall	14"	20"	5"	Short	14.5"
16"	3"	Short	13"	20"	5"	Tall	16"
16"	3"	Tall	14.5"	20"	5"	Short	15"
16"	4"	Short	13.5"	20"	5"	Tall	16.5"
16"	4"	Tall	15"	20"	5"	Short	15.5"
16"	4"	Short	14"	20"	5"	Tall	17"
16"	4"	Tall	15.5"	20"	6"	Short	15.5"
16"	4"	Short	14.5"	20"	6"	Tall	17"
16"	4"	Tall	16"	20"	6"	Short	16"
16"	5"	Short	14.5"	20"	6"	Tall	17.5"
16"	5"	Tall	16"	20"	7"	Short	16.5"
16"	7"	Short	16.5"	20"	7"	Tall	18"
18"	3"	Short	12.5"	22"	4"	Tall	16.5"
18"	3"	Tall	14"				
18"	3"	Short	13"	22"	5"	Tall	16.5"
18"	3"	Tall	14.5"				
18"	4"	Short	13.5"	22"	5"	Tall	17"
18"	4"	Tall	15"				
18"	4"	Short	14"	22"	6"	Tall	17"
18"	4"	Tall	15.5"				
18"	4"	Short	14.5"	22"	6"	Tall	17.5"
18"	4"	Tall	16"				
18"	4"	Short	15"	22"	7"	Tall	18"
18"	4"	Tall	16.5"				
18"	5"	Short	14.5"	24"	6"	Tall	17.5"
18"	5"	Tall	16"				
18"	5"	Short	15"	24"	7"	Tall	18"

Note: All 22" and 24" wheels not with short tilt pivot.

Note: All measurements may vary +/- 3/8"

Note: Measurements may vary +/- 1/2" when choosing foot or hand activated wheel locks.



medical equipment

Form No 17-363CAADL, rev 10/01/2021

NOTE: All specifications and dimensions are approximate
For the most current pricing information visit www.invacare.ca

STROLLER HANDLES

Note: All stroller handle ext. add approximately 6.5" to push handle height.

Note: All stroller/push handle extensions except individual include a black neoprene wrap around pad.

- NXT-30560 Push Handle Extension - Width Adjustable No Charge
- Note: Mounts direct to back post push handle
- NXT-30570 Push Button, Angle - Adjustable Folding Stroller Handle Extension \$121
- Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.
- Note: Mounts direct to back post push handle
- NXT-30575 Push Button, Angle-Adjustable Folding Stroller Handle Extension HD \$193
- Note: Available 10W-22w wheelchair frames.
- Note: Same as above but with metal hinge mechanism.
- Note: Mounts direct to back post push handle.
- NXT-30580 Push Button, Angle-Adjustable, Folding Push Handle Extension \$121
- Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.
- Note: Mounts direct to straight back post tubes.
- NXT-30585 Push Button, Angle-Adjustable Folding Push Handle Extension HD \$193
- Note: Available 10W-22W wheelchair frames.
- Note: Same as above but with metal hinge mechanism.
- Note: Mounts direct to straight back post tubes.
- NXT-30517 Individual Stroller Handle Extensions Removable - Pair \$110



NXT-30560



NXT-30570 / NXT-30575



NXT-30580 / NXT-30585




NXT-30517

LEGRESTS, FOOTRESTS & EXTENSIONS

- NXT-33000 Omit Legrest Hangers & Footrests No Charge
- NXT-30280 75° Pop-Up Swing Away Hangers No Charge
- Note: Not available with high mount bracket.
- NXT-30281 90° Pop-Up Swing Away Hangers No Charge
- Note: Not available with 8" adjustable extensions.
- NXT-30279 90° Pop-Up Swing Away Elevating Legrests (ELRs) \$275
- Note: Not available with all length adjustable extensions.
- Note: Not available with 6.5" heavy duty extensions.
- NXT-30282 75° Latch Style Swing Away Hangers No Charge
- Note: Not available with high mount bracket.
- NXT-30283 90° Latch Style Swing Away Hangers No Charge
- Note: Not available with 8" adjustable extensions.



NXT-17070



NXT-36020

FOOTPLATES & FOOTBOARDS

- NXT-36404 Omit Footrest Plates No Charge
- NXT-17070 Composite Flip-Up No Charge
- NXT-36020 Aluminum Angle Adjustable, Flip-up (14"W frame or narrower) \$173
- NXT-36022 Heel Loops \$44
- NXT-36021 Aluminum Angle Adj. Flip-up w/ Heel Loops (15"W frame or greater) \$194
- o NXT-36023 Omit Heel Loops on NXT-36023 No Charge
- NXT-36024 Calf Strap - Webbing \$44
- NXT-36029 Aluminum Angle Adjustable Footboard Removable \$173



NXT-36040 (MAPS)



NXT-36029

MAPS FOOTPLATE MOUNTING HARDWARE

- NXT-36040 MAPS (Multi Axis Positioning System) - Pair \$330
- NXT-36040-L MAPS (Multi Axis Positioning System) - Left \$165
- NXT-36040-R MAPS (Multi Axis Positioning System) - Right \$165

Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning.

Note: MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.

Note: Must choose aluminum footplates.

Note: Not available on 10" & 11" wide frames

Note: MAPS not available with composite footplates or aluminum footboard.

LEGREST EXTENSIONS (MANDATORY)

- NXT-30285 4" Short Length Adjustable Extensions No Charge
- NXT-30286 6" Medium Length Adjustable Extensions No Charge
- NXT-30287 8" Long Length Adjustable Extensions No Charge

Note: Not available with all 90° hangers.

- NXT-36401 4.5" Heavy Duty ELR Extensions \$60
- NXT-36402 4.5" Heavy Duty Extensions \$60

Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

- NXT-36403 6.5" Heavy Duty Extensions \$71

Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

- NXT-36407 High Mount Brackets \$110

Note: Includes flip-up mounting hardware for future growth adjustment.

Note: Not available with all 75° hangers.

LEGRESTS, FOOTRESTS & EXTENSIONS - SPECIALIZED FOOTREST SYSTEMS

Note: All legrest measurements are from the seat rail to footplate.

Note: Not available with LEAF or SRO seat rail options.

SUPER CONTRACTURE FOOTREST HANGER

- NXT-30291-2 Short Pair - 11" to 15" \$3,852
- NXT-30292-2 Long Pair - 13" to 18" \$3,852

Note: Super Contracture for frames 15" wide and greater.

Note: Super Contracture is both swing-away and removable.

Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5" increments.

Note: Super Contracture has short or long length-adjustable extension tubes with adjustable ankle joints that have 48 positions that adjust in 7.5" increments for plantar / dorsiflex-ion positioning.

Note: Super Contracture includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.

Note: Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference.

Note: Super Contracture position may limit ankle joint/ footplate adjustment range.

3D DYNAMIC FOOTREST HANGER

- NXT-30295-2 Short Pair - 13" to 15" \$4,953
- NXT-30296-2 Long Pair - 15" to 17" \$4,953

Note: 3D Dynamic for frames 15" wide and greater.

Note: 3D Dynamic is both swing-away and removable.

Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.

Note: 3D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel.

Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.

Note: 3D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.

Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.

Footrest / Extension Lengths					
Extension/Hanger	75° Pop	75° Latch	90° Pop	90° Latch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	5"-11"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	12"-16"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	12"-18"
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A
10" Adjustable**	16"-19"	14.5"-18.5"	N/A	N/A	N/A
MAPS Footrest / Extension Lengths					
Extension/Hanger	75° Pop	75° Latch	90° Pop	90° Latch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	11.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	11.5"-17.5"
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A
10" Adjustable**	14.5"-17.5"	14"-18"	N/A	N/A	N/A

Note: Length from top of footrest plate to top of seat rail.

Note: Without extensions bottom hanger hole not accessible for hardware.

Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on he

** 10" Adjustable extension available on PRO-CG ONLY **



NXT-30525

NXT-30530

NXT-30531

Select store location - AADL Numbers

- Calgary - AADL 0000 78000 003
- Edmonton - AADL 0000 78000 001
- Grande Prairie - AADL 0000 665538 001

ARMRESTS & ARMREST PADS

- NXT-38209 Omit Armrests

BACK POST MOUNTED ARMS

Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.

- NXT-38204 10.5" Long Flip-Back No Charge

Note: Full length pad on 13.5" or longer armrests only.

- NXT-38202 13.5" Long Flip-Back No Charge

Note: Full length pad on 13.5" or longer armrests only.

- NXT-38201 10.5" Long Flip-Back Locking \$264

Note: Full length pad on 13.5" or longer armrests only.

- NXT-38200 13.5" Long Flip-Back Locking \$264

Note: Full length pad on 13.5" or longer armrests only.

- NXT-38208 10.5" to 13.5" Long Flip-Back, Telescoping \$259

Note: Full length pad on 13.5" or longer armrests only.

- NXT-38210 10.5" to 13.5" Long Flip-Back, Locking, Telescoping \$286

Note: Full length pad on 13.5" or longer armrests only.

- NXT-38215 17" to 20" Long Flip-Back, Locking, Telescoping \$286

- NXT-38216 Foam Wraps No Charge

- NXT-38213 Desk Length Pads No Charge

- NXT-38214 Full Length Pads No Charge

Note: 8.5" armrests not available with desk and full length armrest pads.

Note: Full length pad on 13.5" or longer armrests only.

Armrest Heights			
Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short	N/A	8" - 12"	8" - 12"
Single Post Height Adj. Tall	N/A	10" - 16"	10" - 16"
Single Post Fxd Hgt Short	N/A	8"-10.5"	8"-10.5"
Single Post Fxd Hgt Tall	N/A	9"-11.5"	9"-11.5"
Dual Post Fixed Height	N/A	9.5"	9.5"
Dual Post Flip Back Height Adj.	N/A	9" - 13"	9" - 13"

Note: All heights are measured from top of seat frame to top of armrest pad.

FRAME MOUNTED ARMRESTS

- NXT-38222 Single Post Fixed Height, Removable - Short \$243

- NXT-38224 Single Post Fixed Height, Removable - Tall \$243

Note: Single Post Fixed Height arms include a removable side guard.

Note: Adjustable in height with tools only. See matrix for height ranges.

- NXT-38217 Single Post Height Adjustable, Removable - Short \$286

- NXT-38203 Single Post Height Adjustable, Removable - Tall \$286

- NXT-38205 Dual Post Fixed Height Removable No Charge

Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

- NXT-38206 Dual Post Flip-Back Height, Adjustable, Removable \$297

Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

- NXT-38213 Desk Length Pads No Charge

- NXT-38214 Full Length Pads No Charge

FRONT WHEELS - CASTERS

- NXT-34333 Omit Caster Wheels No Charge

- NXT-34888 3" Polyurethane \$60

- NXT-34111 3" Polyurethane Light-Up \$60

Note: 3" casters recommended for low seat to floor heights.

- NXT-40000 4" Polyurethane \$71

- NXT-41000 4" Polyurethane Light-Up \$109

- NXT-43000 4" x 1.4" Poly, Aluminum Mag \$99

- NXT-34555 5" x 1" Poly, Plastic Mag No Charge

- NXT-35050 5" Polyurethane Light-Up \$109

- NXT-34222 5" x 1.5" Poly, Plastic Mag \$94

- NXT-42000 5" x 1.4" Poly, Aluminum Mag \$99

- NXT-34666 6" x 1" Poly, Plastic Mag No Charge

- NXT-34999 6" x 1.5" Poly, Plastic Mag \$94

- NXT-44000 6" x 1.4" Poly, Aluminum Mag \$99

- NXT-34444 6" x 2" Poly with Tread \$99

- NXT-34667 6" x 1" Pneumatic Plastic Mag \$83

- NXT-34777 7" x 1" Poly, Plastic Mag \$83

- NXT-34200 7" x 2" Polyurethane \$127

Note: All 7" casters not available with all 12" rear wheels.

Note: All dual post armrests not available with 22" or 24" wheels.

REAR WHEELS, TIRES & ACCESSORIES

Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.

- NXT-31524 Omit Rear Wheels and Axles No Charge

MAG WHEELS

- NXT-31526 12" Low Profile Poly Mag \$39

Note: Not available with any foot/hand activated brake/wheel lock system.

Note: Not available with airless inserts.

Note: All 12" wheels not available with all 7" casters.

- NXT-31523 12" Semi Pneumatic Plastic Mag \$181

Note: Not available with airless inserts.

Note: All 12" wheels not available with all 7" casters.

- NXT-31525 12" Pneumatic Plastic Mag No Charge

- NXT-31019 12" Airless Inserts \$83

Note: All 12" wheels not available with all Handrims.

Note: All 12" wheels not available with all 7" casters.

- NXT-31530 16" Pneumatic Plastic Mag \$44

- NXT-31020 16" Airless Inserts \$83

- NXT-31014 16" Aluminum Handrims \$182

Note: 16" mag wheels not available with plastic coated Handrims.

- NXT-31533 18" Pneumatic Plastic Mag \$44

- NXT-31021 18" Airless Inserts \$83

- NXT-31015 18" Aluminum Handrims \$182

Note: 18" mag wheels not available with plastic coated Handrims.

- NXT-31527 20" Pneumatic Plastic Mag with Aluminum Handrims No Charge

- NXT-31022 20" Airless Inserts \$83

- NXT-31016 20" Plastic Coated Handrim Upgrade \$128

Select store location - AADL Numbers

- Calgary - AADL 0000 78000 003
- Edmonton - AADL 0000 78000 001
- Grande Prairie - AADL 0000 665538 001



medical equipment

Form No 17-363CAADL, rev 10/01/2021

NOTE: All specifications and dimensions are approximate
For the most current pricing information visit www.invacare.ca

REAR WHEELS, TIRES & ACCESSORIES (Continued)

- NXT-31528 22" Pneumatic Plastic Mag with Aluminum Handrims No Charge
- NXT-31023 22" Airless Inserts \$83
- NXT-31017 22" Plastic Coated Handrim Upgrade \$128

Note: All 22" & 24" wheels not available with short tilt bracket.

Note: All 22" & 24" wheels not available with dual post armrests.

- NXT-31529 24" Pneumatic Plastic Mag with Aluminum Handrims No Charge
- NXT-31024 24" Airless Inserts \$83
- NXT-31018 24" Plastic Coated Handrim Upgrade \$128

Note: All 22" & 24" wheels not available with short tilt bracket.

Note: All 22" & 24" wheels not available with dual post armrests.

- NXT-31006 Omit Handrims from 20", 22" or 24" Wheels No Charge

SPOKE WHEELS

- NXT-31534 16" Grey Pneumatic Tire and Silver Spoke Aluminum Wheel No Charge

Note: Rim, hub & spokes are anodized silver.

- NXT-31020 16" Airless Inserts \$83
- NXT-31030 16" Silver Aluminum Handrims \$83
- NXT-31032 16" Black Plastic Coated Handrims \$83

Note: 16" & 18" spoke wheels do not include Handrims.

- NXT-31535 16" Black Pneumatic Tire and Black Spoke Aluminum Wheel No Charge

Note: Rim, hub & spokes are anodized black.

- NXT-31020 16" Airless Inserts \$83
- Note: 16" & 18" spoke wheels do not include Handrims.*
- NXT-31031 16" Black Aluminum Handrims \$83
- NXT-31032 16" Black Plastic Coated Handrims \$83

Note: 16" & 18" spoke wheels do not include Handrims.

- NXT-31537 18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel No Charge

Note: Rim, hub & spokes are anodized silver.

- NXT-31021 18" Airless Inserts \$83
- NXT-31034 18" Silver Aluminum Handrims \$83
- NXT-31036 18" Black Plastic Coated Handrims \$83

Note: 16" & 18" spoke wheels do not include Handrims.

- NXT-31538 18" Black Pneumatic Tire & Black Spoke Aluminum Wheel No Charge

Note: Rim, hub & spokes are anodized black.

- NXT-31021 18" Airless Inserts \$83
- NXT-31035 18" Black Aluminum Handrims \$83
- NXT-31036 18" Black Plastic Coated Handrims \$83

Note: 16" & 18" spoke wheels do not include Handrims.

- NXT-31540 20" Pneumatic Spoke Wheel with Aluminum Handrims \$220
- NXT-31022 20" Airless Inserts \$83

Note: Not available with foot/hand activated wheel locks.

Note: 20", 22" & 24" spoke wheels not available plastic coated Handrims.

- NXT-31541 22" Pneumatic Spoke Wheel with Aluminum Handrims \$220
- NXT-31023 22" Airless Inserts \$83

Note: Not available with foot/hand activated wheel locks.

Note: All 22" & 24" wheels not available with short tilt bracket.

Note: All 22" & 24" wheels not available with dual post armrests.

Note: 20", 22" & 24" spoke wheels not available plastic coated Handrims.

- NXT-31542 24" Pneumatic Spoke Wheel with Aluminum Handrims \$220
- NXT-31024 24" Airless Inserts \$83

Note: Not available with foot/hand activated wheel locks.

Note: All 22" & 24" wheels not available with short tilt bracket.

Note: All 22" & 24" wheels not available with dual post armrests.

Note: 20", 22" & 24" spoke wheels not available plastic coated Handrims.

- NXT-31006 Omit Handrims from 20", 22" or 24" Wheels No Charge

REAR WHEELS, TIRES & ACCESSORIES

- NXT-31501 Push to Lock Wheel Locks (Standard with NXT-30277) No Charge
- NXT-31506 Pull to Lock Wheel Locks No Charge
- NXT-31801 Wheel Lock Extension Handles \$66
- NXT-31508 Foot Activated Wheel Lock (Drum) No Charge

Note: Includes manual wheel locks to be used with public transportation.

Note: Not available with 12" poly plastic mag wheels.

- NXT-30275 Hand Activated Brake and Wheel Lock System (Drum) No Charge

Note: Includes manual wheel locks to be used with public transportation.

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Used to control direction, speed, brake and lock wheels from strolling position.

- NXT-30277 Combination Dual Hand Activated Tilt, Brake and Wheel Lock System (Drum) No Charge

Note: Includes manual wheel locks to be used with public transportation.

Note: Not available with 12" poly plastic mag wheels.

Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position.

REAR ANTI-TIPPERS

- NXT-31600 Curb Assist Anti-Tip Assembly - Flip-Up, Adjustable \$105
- Note: Pair includes one anti-tip tube with curb assist feature.*
- NXT-31603 Standard Anti-Tip Assembly - Flip-Up, Adjustable \$105

SEAT PAN OPTIONS

- NXT-39012 Removable Solid Seat Pan - Bolteron \$183
- NXT-39032 Removable Solid Seat Pan - Aluminum \$275

Note: Removable Seat Pan raises the effective seat to floor height by 5/8".

Note: Aluminum Removable Solid Seat Pan adjusts 2" in width by mounting hardware to grow with NXT frame.

Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.

Note: All seat pans include 2" x 2" notches for growth adjustment.

Note: All seat pans made to fit frame dimensions unless specified below.

Seat Width _____" X Seat Depth _____"

Note: Please note desired seat pan size. 4" extra depth maximum.

Select store location - AADL Numbers

- Calgary - AADL 0000 78000 003
- Edmonton - AADL 0000 78000 001
- Grande Prairie - AADL 0000 665538 001



ACCESSORIES

<input type="checkbox"/>	NXT-23207	Medical Necessities Pouch - Black Nylon	\$73
<input type="checkbox"/>	NXT-SFDI-9HP	Black Neoprene Frame Protectors - For Hangers	\$99
<input type="checkbox"/>	NXT-10911-2	Black Neoprene Frame Protectors - For Back Posts	\$88
<input type="checkbox"/>	NXT-23210	Flip-Down Rigidizing Bar	\$231

Note: Mounts to back post tubes to add rigidity to frame.

Note: Does NOT adjust in width to grow with frame.

<input type="checkbox"/>	NXT-10900	Wheelchair Backpack	\$60
<input type="checkbox"/>	NXT-10905	Wheelchair Canopy w/ Window and Storage Pocket - Flip-Back	\$138

Note: Canopy assembly includes removable mounting hardware.

ANTERIOR SUPPORT SURFACES - TRAYS

<input type="checkbox"/>	NXT-08-2000	3/8" Lexan Tray - Channel Bracket Hardware	\$292
<input type="checkbox"/>	NXT-08-6000	3/8" Lexan Tray - Adjustable Knob Hardware & Ribs	\$236
<input type="checkbox"/>	NXT-08-13000	1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs	\$197
<input type="checkbox"/>	NXT-08-16000	1/4" Lexan Tray - Hook & Loop Straps	\$114
<input checked="" type="checkbox"/>	NXT-10917	Wheelchair Tool Kit	Standard

SEATING

Note: All seating will ship separately of the chair.

Note: Not all seating options are listed. See the full Matrix seating order form for more options.

Note: Prices listed apply to standard sizes only unless otherwise noted.

BACK SUPPORTS

MATRIX MINI - PEDIATRIC POSTURE BACK

<input type="checkbox"/>	MINIXXXX-S	Black Startex Fabric Cover	\$845
<input type="checkbox"/>	MINIXXXX-K	Kid*a*bra Meshtex Fabric Cover	\$845
<input type="checkbox"/>	MINIXXXX-B	Blue Blizzard Startex Fabric Cover	\$845
<input type="checkbox"/>	MINIXXXX-P	Pink Lollipop Startex Fabric Cover	\$845

ELITE TR - TILT & RECLINE BACK WITH SUPERSOFT FOAM

Provides an additional 1" layer of high resiliency foam for added comfort and immersion.

<input type="checkbox"/>	PBEXXXX-TR	Elite TR Back	\$866
<input type="checkbox"/>	PBEWSCXXX-TR	Wipeable Slipcover Add-On for Elite TR Back (Stock sizes only)	\$150

ELITE - POSITIONING BACK

<input type="checkbox"/>	PBEXXXX	PB Elite Back	\$766
<input type="checkbox"/>	PBEWSCXXX	Wipeable Slipcover Add-On for Elite Back (Stock sizes only)	\$150

CUSHIONS

KID*A*BRA VI - PEDIATRIC SKIN PROTECTION AND POSITIONING FOAM CUSHION

<input type="checkbox"/>	KAXXXX-VI	Kid*a*bra Meshtex Fabric Cover	\$545
--------------------------	-----------	--------------------------------	-------

MINI PS - PEDIATRIC POSTURE SEAT FOAM CUSHION

<input type="checkbox"/>	PSXXXX	Reversible Black Startex Cover	\$569
<input type="checkbox"/>	PSXXXX-K	Kid*a*bra Startex Fabric Cover (Non-reversible)	\$569
<input type="checkbox"/>	PSXXXX-B	Blue Blizzard Startex Fabric Cover (Non-reversible)	\$569
<input type="checkbox"/>	PSXXXX-P	Pink Lollipop Startex Fabric Cover (Non-reversible)	\$569

LIBRA - SKIN PROTECTION AND POSITIONING FOAM AND FLUID CUSHION

<input type="checkbox"/>	LCXXXX	Startex Cover, Fabric Side Out <= 20W x 20D	\$745
<input type="checkbox"/>	LCXXXX-IC	Infection Control Cover <= to 20W x 20D	\$798

ACCESSORIES

BASIC HEADREST

<input type="checkbox"/>	5500071	Comfort Headrest Assembly w/ Multi Axis Hardware	\$375
--------------------------	---------	--	-------

ELAN HEADREST PADS - HARDWARE SOLD SEPARATELY

<input type="checkbox"/>	ESP6	Standard Pad - Startex - tailored, removable cover. 6" x 3.5"	\$310
<input type="checkbox"/>	ESP10	Standard Pad - Startex - tailored, removable cover. 10" x 5"	\$310
<input type="checkbox"/>	ESP14	Standard Pad - Startex - tailored, removable cover. 14" x 5"	\$310
<input type="checkbox"/>	ESP6-IC	Standard Pad - Infection Control - tailored, removable cover. 6" x 3.5"	\$335
<input type="checkbox"/>	ESP10-IC	Standard Pad - Infection Control - tailored, removable cover. 10" x 5"	\$335
<input type="checkbox"/>	ESP14-IC	Standard Pad - Infection Control - tailored, removable cover. 14" x 5"	\$335
<input type="checkbox"/>	EOP9	Occipital Pad - Startex - tailored, removable cover. 9" x 7"	\$341
<input type="checkbox"/>	EOP12	Occipital Pad - Startex - tailored, removable cover. 12" x 8"	\$341
<input type="checkbox"/>	EOP9-IC	Occipital Pad - Infection Control - tailored, removable cover. 9" x 7"	\$366
<input type="checkbox"/>	EOP12-IC	Occipital Pad - Infection Control - tailored, removable cover. 12" x 8"	\$366
<input type="checkbox"/>	E4POINT	4 Point Positioning Pad - Startex - tailored, removable cover	\$448
<input type="checkbox"/>	E4POINT-IC	4 Point Positioning Pad - Infection Control - tailored, removable cover	\$473

ELAN HEADREST HARDWARE - PADS SOLD SEPARATELY

<input type="checkbox"/>	MEHW	Standard Elan Hardware - 15" Vertical Post	\$266
<input type="checkbox"/>	MEMH	Mini Elan Hardware - 11" Vertical Post	\$254



Kid*a*bra Meshtex Fabric



Blue Blizzard Startex Fabric



Pink Lollipop Startex Fabric

NOTES:



Select store location - AADL Numbers

- Calgary - AADL 0000 78000 003
- Edmonton - AADL 0000 78000 001
- Grande Prairie - AADL 0000 665538 001

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)

Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

NOTE: All specifications and dimensions are approximate