

Select store location - AADL Numbers

- ☐ Calgary AADL 00000 78000 003
- ☐ Edmonton AADL 00000 78000 001 ☐ Grande Prairie - AADL 0000 665538 001



Freedom® NXT™ Folding Tilt-in-Space Wheelchair **AADL Customers Only AADL # W684** Effective October 1st, 2021

| | | PP-1 | |
|-----------|---------------|-----------|--|
| NAME | | PHONE NO_ | |
| ACCOUNT # | CUSTOMER P.O. | DATE | |
| ADDRESS | | | |

*This chair will ship directly to the address specified from Freedom Designs.

FRAME TYPE

■ NXT-30301 NXT Pediatric Folding Tilt-in-Space Wheelchair w/ Transport Package \$3,650

Note: Includes Freedom Future Fit Program

TRANSPORT TIE DOWN

■ NIXT_00000_1 Transportation Package - For NXT Pediatric Standard

SPECIALIZED SEAT RAIL OPTIONS

Note: Not required for standard swingaway or pop-up footrests.

Note: Not available with Contracture & 3D Dynamic footrest systems.

| LEAF: LOWER EXTR | REMITY ACCOMODATION FRAME | |
|----------------------------|---|-------|
| □ NXT-30324 | NXT LEAF Seat Frame - Pair | \$770 |
| □ NXT-30324-L | NXT LEAF Seat Frame - Left | \$385 |
| ■ NXT-30324-R | NXT LEAF Seat Frame - Right | \$385 |
| Note: LEAF front frame and | gle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. | |
| | | |

Note: LEAF seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.

Note: No seating system is provided with NXT frames. A seating system must be purchased separately.

SRO: SEAT RAIL ONLY

| ■ NXT-30325 | NXT SRO Seat Frame- Pair | \$330 |
|---------------|---------------------------|-------|
| □ NXT-30325-L | NXT SRO Seat Frame- Left | \$165 |
| ■ NXT-30325-R | NXT SRO Seat Frame- Right | \$165 |

Note: SRO is manufactured without standard footrest hanger attachments.

Note: SRO is used when customer attaches leg/footrest(s) directly to the underside of the solid seat.

TILTING OPTIONS

FRAM

| □ NXT-30272 | Dual Lever Hand Activated Tilt | Standard |
|-------------|---|-----------|
| □ NXT-30273 | Foot Activated Tilt | No Charge |
| □ NXT-30277 | Combination Dual Hand Activated Tilt, Brake and Wheel Lock System | No Charge |

Note: Comes standard with manual push to lock wheel locks, NXT-31501. Pull to lock optional.

Note: Not available with 12" Low Profile Poly Mag Wheels

| 12004 | Split Collar Mechanism | Standard |
|-------------------|---|-----------|
| ME FINISHES | | |
| □ SUGARCOOKIE | Sugar Cookie | No Charge |
| □ NEONY | Neon Yellow | No Charge |
| ☐ YELLOWY | Yellow | No Charge |
| □ BRIGHTP | Bright Pink | No Charge |
| ☐ PINKP | Pink Parfait | No Charge |
| □ SUGARP | Sugar Plum | No Charge |
| ☐ CANDYV | Candy Violet | No Charge |
| ☐ CANDYR | Candy Raspberry | No Charge |
| □ PACIFICB | Pacific Blue | No Charge |
| ☐ BLUEB | Blue Blast | No Charge |
| □ ELECTRICB | Electric Blue | No Charge |
| ☐ GALAXYB | Galaxy Blue | No Charge |
| ☐ PASSIONP | Passion Purple | No Charge |
| □ OCTANEO | Octane Orange | No Charge |
| ☐ REDS | Red Sunset | No Charge |
| ☐ REALLYR | Really Red | No Charge |
| □ EXTREMEG | Extreme Green | No Charge |
| ☐ FIERYR | Fiery Red | No Charge |
| □ SPARKLEG | Sparkle Green | No Charge |
| ☐ TWILIGHTT | Twilight | No Charge |
| □ MYSTICT | Mystic Teal | No Charge |
| □ STORMYG | Stormy Gray | No Charge |
| ☐ BLACKB | Black | No Charge |
| WIDTH RANGE SETTI | NGS Note: Frame widths measured from outside of seat rails. | |
| □ NXT-30250 | 10" / 11" / 12" Width Range | No Charge |

SEAT

| | BLACKB | Black | No Charge |
|-----|---------------|--|-----------|
| WID | TH RANGE SETT | INGS Note: Frame widths measured from outside of seat rails. | Ü |
| | NXT-30250 | 10" / 11" / 12" Width Range | No Charge |
| 0 | NXT-W10 | Set at 10" Wide | |
| 0 | NXT-W11 | Set at 11" Wide | |
| 0 | NXT-W12 | Set at 12" Wide | |
| | NXT-30210 | 11" / 12" / 13" Width Range | No Charge |
| 0 | NXT-W11 | Set at 11" Wide | |
| 0 | NXT-W12 | Set at 12" Wide | |
| 0 | NXT-W13 | Set at 13" Wide | |
| | NXT-30220 | 12" / 13" / 14" Width Range | No Charge |
| 0 | NXT-W12 | Set at 12" Wide | |
| 0 | NXT-W13 | Set at 13" Wide | |
| 0 | NXT-W14 | Set at 14" Wide | |
| | NXT-30230 | 13" / 14" / 15" Width Range | No Charge |
| 0 | NXT-W13 | Set at 13" Wide | |
| 0 | NXT-W14 | Set at 14" Wide | |
| 0 | NXT-W15 | Set at 15" Wide | |
| | NXT-30240 | 14" / 15" / 16" Width Range | No Charge |
| 0 | NXT-W14 | Set at 14" Wide | |
| 0 | NXT-W15 | Set at 15" Wide | |
| 0 | NXT-W16 | Set at 16" Wide | |
| | NXT-30290 | 15" / 16" / 17" Width Range | No Charge |
| 0 | NXT-W15 | Set at 15" Wide | |
| 0 | NXT-W16 | Set at 16" Wide | |
| 0 | NXT-W17 | Set at 17" Wide | |



Weight Capacity: 225 lbs Limited 5 Year Warranty on Frame and Axles

| NXT OVERALL CHAIR WIDTH | | | | | |
|-------------------------|------|----------------------|---------------|--|--|
| FRAME WIDTH + | | | | | |
| W | HEEL | STD WHEEL LOCK | DRUM BRAKE | | |
| | 12" | +7.5" | +8.5" | | |
| | 16" | +7.5" | +9.5" | | |
| MAG | 18" | +7.5" | +9.5" | | |
| IVIAG | 20" | +7.5" | +8.5" | | |
| | 22" | +7.5" | +8.5" | | |
| | 24" | +7.5" | +8.5" | | |
| | 16" | +7.5" | +8.5" | | |
| | 18" | +7.5" | +8.5" | | |
| SPOKE | 20" | +7.5" | +8.5" | | |
| | 22" | +7.5" | +8.5" | | |
| | 24" | +7.5" | +8.5" | | |

| FREEDOM NX1 | Г™ Ас | cessories and Option | าร | | | | | <u>N</u> 0 | <u>OTES</u> | | | |
|---------------------------------------|---------------|---|---|------------------------|------------|--------------|---------------|---------------|--------------|--------------|-------------|------------|
| SEAT WIDTH RANGE | | ` ' | | | | | | | | | | |
| Note: Frame widths measured NXT-30260 | 0 | 16" / 17" / 18" Width Range | | No Charge | | | | | | | | |
| 0 NXT-V 0 NXT-V | | Set at 16" Wide Set at 17" Wide | | | | | | | | | | |
| o NXT-V | W18 | Set at 18" Wide | | | | | | | | | | |
| □ NXT-3026 ⁻ | 61 | 17" / 18" Width Range | | No Charge | | | | | | | | |
| o NXT-V | W17 | Set at 17" Wide | | ū | | | | | | | | |
| o NXT-V | W18 | Set at 18" Wide | | | | | | | | | | |
| ☐ NXT-30262 | 52 | 18" - Not width adjustable | | No Charge | | | | | | | | |
| o NXT-V | W18 | Set at 18" Wide | | | | | | | | | | |
| SEAT DEPTH RANGI | E SETTII | NGS | | | | | | | | | | |
| | | | om front of back post to end of the seat frame. | For back post lo | ocation ad | d seating / | hardware | thickness | (typically 2 | ") to client | 's actual s | eat depth. |
| | | hardware thickness of 2" = back po | | | | | | | | | | |
| | | | ted positions. Adjusting seat rail location will described by | | | imitations | and possi | ble instabil | ity. | | | |
| □ NXT-3051 | | SHORT FRAME & SHORT TI | LT PIVOT BRACKET | No Charge | | NX | T Wheel | / Caster | Seat to F | loor Heig | hts | |
| o NXT-30 | | Set at 12" Depth | | No Charge | Wheel | Caster | Pivot | Height | Wheel | Caster | Pivot | Height |
| 0 NXT-30 | | Set at 13" Depth | | No Charge | 12" | 3" | Short | 12.5" | 18" | 5" | Tall | 16.5" |
| 0 NXT-30 | | Set at 14" Depth | | No Charge | 12" | 3" | Tall | 14" | 18" | 5" | Short | 15.5" |
| 0 NXT-30 | | Set at 15" Depth | | No Charge | 12" | 3" | Short | 13" | 18" | 5" | Tall | 17" |
| 0 NXT-30 | | Set at 16" Depth | | No Charge | 12" | 3" | Tall | 14.5" | 18" | 6" | Short | 15.5" |
| 0 NXT-30 | | 17" - Must select 2" extension | | No Charge | 12" | 4" | Short | 13.5" | 18" | 6" | Tall | 17" |
| o NXT-30 | | 18" - Must select 2" extension | | No Charge | | | | | | | | |
| ☐ NXT-3052 | | 2" Extension - Pair | | \$33 | 12" | 4" | Tall | 15" | 20" | 4" | Short | 14" |
| | | ight range is 12.5 to 16.5". | | | 12" | 5" | Short | 15.5" | 20" | 4" | Tall | 15.5" |
| | | 9" and 20" not available with short fra | | | 12" | 5" | Tall | 17" | 20" | 4" | Short | 14.5" |
| | | available with short frame and short tilt | | N. O. | 12" | 6" | Short | 16" | 20" | 4" | Tall | 16" |
| □ NXT-30512 | | MEDIUM FRAME & SHORT T | ILI PIVOI BRACKEI | No Charge | 12" | 6" | Tall | 17.5" | 20" | 4" | Short | 15" |
| 0 NXT-30 | | 14" | | No Charge | 16" | 3" | Short | 12.5" | 20" | 4" | Tall | 16.5" |
| 0 NXT-30 | | 15" | | No Charge | 16" | 3" | Tall | 14" | 20" | 5" | Short | 14.5" |
| 0 NXT-30 | | 16" | | No Charge | 16" | 3" | Short | 13" | 20" | 5" | Tall | 16" |
| 0 NXT-30 | | 17" | | No Charge | 16" | 3" | Tall | 14.5" | 20" | 5" | Short | 15" |
| 0 NXT-30 | | 18" | | No Charge | 16" | 4" | Short | 13.5" | 20" | 5" | Tall | 16.5" |
| 0 NXT-30 | | 19" - Must select 2" extension | | No Charge | 16" | 4" | Tall | 15" | 20" | 5" | Short | 15.5" |
| 0 NXT-30 | | 20" - Must select 2" extension | | \$625 | 16" | 4" | Short | 14" | 20" | 5" | Tall | 17" |
| 0 NXT-30 | | 21" - Must select 3" extension | | \$625 | 16" | 4" | Tall | 15.5" | 20" | 6" | Short | 15.5" |
| □ NXT-30525 | | 2" Extension - Pair | | \$33 | 16" | 4" | Short | 14.5" | 20" | 6" | Tall | 17" |
| □ NXT-30530 | | 3" Extension - Pair | | \$33 | 16" | 4" | Tall | 16" | 20" | 6" | Short | 16" |
| | | ight range is 12.5 to 16.5". | | | 16" | 5" | Short | 14.5" | 20" | 6" | Tall | 17.5" |
| | | available with medium frame and short | | No Charge | 16" | 5" | Tall | 16" | 20" | 7" | Short | 16.5" |
| □ NXT-30540 | | MEDIUM FRAME & TALL TIL | I PIVOI BRACKET | No Charge | 16" | 7" | Short | 16.5" | 20" | 7" | Tall | 18" |
| 0 NXT-305 | | 14" | | No Charge | 16" | 7" | Tall | 18" | 20 | | Iall | 10 |
| 0 NXT-305 | | 15" | | No Charge | 18" | 3" | Short | 12.5" | 22" | 4" | Tall | 16.5" |
| 0 NXT-305 | | 16" 17" | | No Charge | | 3" | | | | 4 | rall | 10.5 |
| 0 NXT-305 | | 18" | | No Charge | | | Tall | 14" | 20" | | T-II | 40.5" |
| 0 NXT-305 0 NXT-305 | | 19" - Must select 2" extension | | No Charge No Charge | 18" 18" | 3" 3" | Short Tall | 13" | 22" | 5" | Tall | 16.5" |
| 0 NXT-305 | | 20" - Must select 2" extension | | \$625 | | | | 14.5" | 00" | 5 " | T. II | 470 |
| 0 NXT-305 | | 21" - Must select 3" extension | | \$625 | 18" | 4" | Short | 13.5" | 22" | 5" | Tall | 17" |
| 0 NXT-305 | | 22" - Must select 4" extension | | \$625 | 18" | 4" | Tall | 15" | 001 | 0" | - " | 470 |
| □ NXT-30525 | | 2" Extension - Pair | | \$33 | 18" | 4" | Short | 14" | 22" | 6" | Tall | 17" |
| □ NXT-30530 | | 3" Extension - Pair | | \$33 | 18" | 4" | Tall | 15.5" | | | | |
| □ NXT-3053 | | 4" Extension - Pair | | \$33 | 18" | 4" | Short | 14.5" | 22" | 6" | Tall | 17.5" |
| | | | | φοο | 18" | 4" | Tall | 16" | | | | |
| □ NXT-30276 | | ight range is 14 to 18". 5° Anterior Tilt Seat Frame | | No Charge | 18" | 4" | Short | 15" | 22" | 7" | Tall | 18" |
| | | 10°, includes split collars for adjustmen | t | No Charge | 18" | 4" | Tall | 16.5" | | | | |
| | | nd Long frames only. | | | 18" | 5" | Short | 14.5" | 24" | 6" | Tall | 17.5" |
| | | • | red from front of back post to end of the seat fram | | 18" | 5" | Tall | 16" | | | | |
| | | | • | ie. | 18" | 5" | Short | 15" | 24" | 7" | Tall | 18" |
| • | | eating / hard- ware thickness (typically | | | | 2" and 24" w | | | pivot. | | | |
| ⊏xampie: seat de | eptn of 12" · | + back / hardware thickness of 2" = ba | ck post location set at 14" | | | neasurement | | | | | | |
| Seat to Flo | oor Heigh | t Requested | " *REQUIRED* | | Note: Mea | surements n | nay vary +/- | - 1/2" when o | hoosing foot | or hand act | ivated whee | l locks. |
| PUSH HANDLES | oor rieigii | | | | | | | | | | | |
| □ NXT-30518 | 8 | Standard 90° Bend Push Hand | dles and Back Post Tubes | No Charge | | | | | | | | |
| □ NXT-30519 | | Ergo Bend Push Handles and | | No Charge | | | | | | | | |
| | | e extensions not recommended with E | | 590 | | | | | | | | |
| □ NXT-3053 | | Seating Dynamics Rocker Bac | | \$2 168 | | | | | | | | |

☐ NXT-30532 Seating Dynamics Rocker Back Note: Angle adjustable from 72° to 108° in 6° increments.

Note: Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.

Note: Locks into a non-dynamic position when necessary. Note: Does not allow the back post tubes to fold forward.

BACK POST HEIGHT

Note: Back Posts are measured from top of seat rail to top of handle bend.

NXT-30515 Short Push Handle - 20" high

NXT-30514 Medium Push Handle - 23" high

NXT-30516 Tall Push Handle - 26" high No Charge No Charge No Charge



Select store location - AADL Numbers □ Calgary - AADL 00000 78000 003

| | Edmonton - AADL 00000 78000 001 |
|--------|-------------------------------------|
| \Box | Grande Prairie - AADI 0000 665538 0 |

STROLLER HANDLES Note: All stroller handle ext. add approximately 6.5" to push handle height. Note: All stroller/push handle extensions except Individual include a black ned □ NXT-30560 Push Handle Extension - Width Adjustable No Charge □ NXT-30570 Push Button, Angle - Adjustable Folding Stroller Handle Extension \$121 Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below e: Mounts direct to back post push handle \$193 Push Button, Angle-Adjustable Folding Stroller Handle Extension HD Note: Available 10W-22w wheelchair frames. Note: Same as above but with metal hinge mechanism. e: Mounts direct to back post push handi Push Button, Angle-Adjustable, Folding Push Handle Extension \$121 □ NXT-30580 Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below Note: Mounts direct to straight back post tubes. \$193 □ NXT-30585 Push Button, Angle-Adjustable Folding Push Handle Extension HD Note: Available 10W-22W wheelchair frames. Note: Same as above but with metal hinge mechanism Note: Mounts direct to straight back post tubes.

| I ECDESTS | FOOTRESTS & EXTENSIONS | |
|-----------|------------------------|--|

| | | * |
|--------------------------------|---|-----------|
| RESTS, FOOTRESTS & E | XTENSIONS | |
| □ NXT-33000 (| Omit Legrest Hangers & Footrests | No Charge |
| □ NXT-30280 7 | 75° Pop-Up Swing Away Hangers | No Charge |
| Note: Not available with high | mount bracket. | |
| □ NXT-30281 9 | 90° Pop-Up Swing Away Hangers | No Charge |
| Note: Not available with 8" ad | ljustable extensions. | |
| □ NXT-30279 9 | 90° Pop-Up Swing Away Elevating Legrests (ELRs) | \$275 |

Individual Stroller Handle Extensions Removable - Pair

Note: Not available with all length adjustable extensions. Note: Not available with 6.5" heavy duty extension

No Charge □ NXT-30282 75° Latch Style Swing Away Hangers e: Not available with high mount bracket □ NXT-30283 90° Latch Style Swing Away Hangers No Charge

Note: Not available with 8" adjustable extensions.

FOOTPLATES & FOOTBOARDS

□ NXT-30517

| □ NXT-36404 | Omit Footrest Plates | No Charge |
|------------------|---|-----------|
| □ NXT-17070 | Composite Flip-Up | No Charge |
| □ NXT-36020 | Aluminum Angle Adjustable, Flip-up (14"W frame or narrower) | \$173 |
| □ NXT-36022 | Heel Loops | \$44 |
| □ NXT-36021 | Aluminum Angle Adj. Flip-up w/ Heel Loops (15"W frame or greater) | \$194 |
| o NXT-36023 | Omit Heel Loops on NXT-36023 | No Charge |
| □ NXT-36024 | Calf Strap - Webbing | \$44 |
| □ NXT-36029 | Aluminum Angle Adjustable Footboard Removable | \$173 |
| FOOTEL ATE MOUNT | TING HADDWADE | |

MAPS FOOTPLATE MOUNTING HARDWARE

| □ NXT-36040 MAPS (Multi Axis Positioning System) - Pair | \$330 |
|--|-------|
| □ NXT-36040-L MAPS (Multi Axis Positioning System) - Left | \$165 |
| □ NXT-36040-R MAPS (Multi Axis Positioning System) - Right | \$165 |
| Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. | ψισσ |

Note: MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning

Note: Must choose aluminum footplates. Note: Not available on 10" & 11" wide frames

Note: MAPS not available with composite footplates or aluminum footboard.

LEGREST EXTENSIONS (MANDATORY) □ NXT-30285 4" Short Length Adjustable Extensions

| ☐ NXT-30286 | 6" Medium Length Adjustable Extensions | No Charge |
|------------------------------|---|-----------|
| ☐ NXT-30287 | 8" Long Length Adjustable Extensions | No Charge |
| Note: Not available with all | 90° hangers. | |
| □ NXT-36401 | 4.5" Heavy Duty ELR Extensions | \$60 |
| □ NXT-36402 | 4.5" Heavy Duty Extensions | \$60 |
| Note: Heavy duty extension | s do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tul | bes. |
| □ NXT-36403 | 6.5" Heavy Duty Extensions | \$71 |
| Note: Heavy duty extension | is do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tul | bes. |

High Mount Brackets Note: Includes flip-up mounting hardware for future growth adjustment.

Note: Not available with all 75° hange

□ NXT-36407

LEGRESTS, FOOTRESTS & EXTENSIONS - SPECIALIZED FOOTREST SYSTEMS

Note: All legrest measurements are from the seat rail to footplate.

Note: Not available with LEAF or SRO seat rail options. SUPER CONTRACTURE FOOTREST HANGER

| OU EN CONTRACTORE FOOTREST HANGER | | | | |
|-----------------------------------|-------------------------|---------|--|--|
| □ NXT-30291-2 | Short Pair - 11" to 15" | \$3,852 | | |
| □ NXT-30292-2 | Long Pair - 13" to 18" | \$3,852 | | |

Note: Super Contracture for frames 15" wide and greater Note: Super Contracture is both swing-away and removable

Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments.

Note: Super Contracture has short or long length-adjustable extension tubes with adjustable ankle joints that have 48 positions that adjust in 7.5° increments for plantar / dorsiflex-ion positioning

Note: Super Contracture includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.

Note: Frame and/or front caster choice may limit Super Contrac- ture adjustment range and/or cause footplate/front caster interference

Note: Super Contracture position may limit ankle joint/ footplate adjustment range.

3D DYNAMIC FOOTREST HANGER

| □ NXT-30295-2 | Short Pair - 13" to 15" | \$4,953 |
|---------------|-------------------------|---------|
| □ NXT-30296-2 | Long Pair - 15" to 17" | \$4,953 |

Note: 3D Dynamic for frames 15" wide and greater

Note: 3D Dynamic is both swing-away and removable.

Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.

Note: 3D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel.

Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning

Note: 3D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.

Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference



NXT-30560

\$110

No Charge

\$110



NXT-30570 / NXT-30575



NXT-30580 / NXT-30585



NXT-30517



NXT-17070



NXT-36020



NXT-36040 (MAPS)



NXT-36029

| Footrest / Extension Lengths | | | | | |
|------------------------------------|-----------|-------------|-----------|-----------|------------|
| Extension/Hanger 75° Pop 75° Latch | | | 90° Pop | 90 °Latch | ELR (90°) |
| Hanger Only | 6"-8" | 5.5"-8.5" | 4"-8" | 4"-8" | 5"-11" |
| High Mount | N/A | N/A | 1.5"-3" | N/A | N/A |
| 4.5" H D Ext. | 9"-12" | 9.5"-12.5" | 9"-12" | 9"-12" | 12"-16" |
| 6.5" H D Ext. | 9"-15" | 9.5"-15.5" | 9"-15" | 9"-15" | 12"-18" |
| 4" Adjustable | 10"-13" | 9.5"-12.5" | 11"-13" | 10"-13" | N/A |
| 6" Adjustable | 12"-15" | 10.5"-14.5" | 13"-15" | 12"-15" | N/A |
| 8" Adjustable | 14"-17" | 12.5"-16.5" | N/A | N/A | N/A |
| 10" Adjustable** | 16"-19" | 14.5"-18.5" | N/A | N/A | N/A |
| MAPS Footrest / Extension Lengths | | | | | |
| Extension/Hanger | 75° Pop | 75° Latch | 90° Pop | 90 °Latch | ELR (90°) |
| Hanger Only | 5.5"-7.5" | 5"-8" | 5.5"-7.5" | 4.5"-7.5" | 4.5"- 8.5" |
| | | | | | |

| <u>.</u> | | | | | |
|------------------|-------------|-----------|-------------|------------|-------------|
| Extension/Hanger | 75° Pop | 75° Latch | 90° Pop | 90 °Latch | ELR (90°) |
| Hanger Only | 5.5"-7.5" | 5"-8" | 5.5"-7.5" | 4.5"-7.5" | 4.5"- 8.5" |
| High Mount | N/A | N/A | N/A | N/A | N/A |
| 4.5" H D Ext. | 8.5"-11.5" | 9"-12" | 8.5"-11.5" | 8.5"-11.5" | 11.5"-15.5" |
| 6.5" H D Ext. | 8.5"-13.5" | 9"-14" | 8.5"-13.5" | 8.5"-13.5" | 11.5"-17.5 |
| 4" Adjustable | 8.5"-11.5" | 9"-12" | 9.5"-11.5" | 8.5"-11.5" | N/A |
| 6" Adjustable | 10.5"-13.5" | 10"-14" | 11.5"-13.5" | 8.5"-13.5" | N/A |
| 8" Adjustable | 12.5"-15.5" | 12"-16" | N/A | N/A | N/A |
| 10" Adjustable** | 14.5"-17.5" | 14"-18" | N/A | N/A | N/A |
| | | | | | |


Note: Length from top of footrest plate to top of seat rail

Note: Without extensions bottom hanger hole not accessible for hardware

Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on he ** 10" Adjustable extension available on PRO-CG ONLY **

NXT-30525 NXT-30530

NXT-30531



- Edmonton AADL 00000 78000 001
 - Grande Prairie AADL 0000 665538 001



ARMRESTS & ARMREST PADS

Omit Armrests □ NXT-38209

| RACK | PO: | STN | MOUNTED | ARMS |
|------|-----|-----|---------|------|

| DAGITI COI MICCITI | LD ARMO | |
|-----------------------------|--|-----------|
| | ngth is measured from the front of the backpost tube to the end of the armrest with foam | |
| □ NXT-38204 | 10.5" Long Flip-Back | No Charge |
| • • | 2.5" or longer armrests only. | |
| □ NXT-38202 | 13.5" Long Flip-Back | No Charge |
| Note: Full length pad on 13 | 3.5" or longer armrests only. | |
| □ NXT-38201 | 10.5" Long Flip-Back Locking | \$264 |
| Note: Full length pad on 13 | 2.5" or longer armrests only. | |
| ☐ NXT-38200 | 13.5" Long Flip-Back Locking | \$264 |
| Note: Full length pad on 13 | 3.5" or longer armrests only. | |
| □ NXT-38208 | 10.5" to 13.5" Long Flip-Back, Telescoping | \$259 |
| Note: Full length pad on 13 | 3.5" or longer armrests only. | |
| □ NXT-38210 | 10.5" to 13.5" Long Flip-Back, Locking, Telescoping | \$286 |
| Note: Full length pad on 13 | 3.5" or longer armrests only. | |
| □ NXT-38215 | 17" to 20" Long Flip-Back, Locking, Telescoping | \$286 |
| □ NXT-38216 | Foam Wraps | No Charge |
| □ NXT-38213 | Desk Length Pads | No Charge |
| □ NXT-38214 | Full Length Pads | No Charge |
| Note: 8.5" armrests not ava | ailable with desk and full length armrest pads. | |
| AL C = 0.1 0 1 40 | | |

| Armrest Heights | | | | |
|--|--------------|--------------|--------------|--|
| Armrest and Pad | Foam | Desk | Full | |
| Flip Back | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" | |
| Flip Back Locking | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" | |
| Telescoping Flip Back | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" | |
| Flip Back Locking Telescoping | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" | |
| Single Post Height Adj. Short | N/A | 8" - 12" | 8" - 12" | |
| Single Post Height Adj. Tall | N/A | 10" - 16" | 10" - 16" | |
| Single Post Fxd Hgt Short | N/A | 8"-10.5" | 8"-10.5" | |
| Single Post Fxd Hgt Tall | N/A | 9"-11.5" | 9"-11.5" | |
| Dual Post Fixed Height | N/A | 9.5" | 9.5" | |
| Dual Post Flip Back Height Adj. | N/A | 9" - 13" | 9" - 13" | |
| Note: All heights are measured from top of seat frame to top of armrest pad. | | | | |

FRAME MOUNTED ARMRESTS

Note: Full length pad on 13.5" or longer armrests only.

| □ NXT-38222 | Single Post Fixed Height, Removable - Short | \$243 |
|------------------|---|-------|
| □ NXT-38224 | Single Post Fixed Height, Removable - Tall | \$243 |
| N (0: (D (E:) | | |

Note: Single Post Fixed Height arms include a removable side guard. Note: Adjustable in height with tools only. See matrix for height ranges

Single Post Height Adjustable, Removable - Short ☐ NXT-38217 \$286 Single Post Height Adjustable, 1561.

Dual Post Fixed Height Removable □ NXT-38203 Single Post Height Adjustable, Removable - Tall \$286 □ NXT-38205 No Charge

Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

□ NXT-38206 Dual Post Flip-Back Height, Adjustable, Removable \$297 Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

Desk Length Pads □ NXT-38213 No Charge □ NXT-38214 Full Length Pads No Charge

FRONT WHEELS - CASTERS

| □ NXT-34333 | Omit Caster Wheels | No Charge | | |
|--|--------------------------------------|-----------|--|--|
| □ NXT-34888 | 3" Polyurethane | \$60 | | |
| □ NXT-34111 | 3" Polyurethane Light-Up | \$60 | | |
| Note: 3" casters recomm | ended for low seat to floor heights. | | | |
| □ NXT-40000 | 4" Polyurethane | \$71 | | |
| ■ NXT-41000 | 4" Polyurethane Light-Up | \$109 | | |
| □ NXT-43000 | 4" x 1.4" Poly, Aluminum Mag | \$99 | | |
| ■ NXT-34555 | 5" x 1" Poly, Plastic Mag | No Charge | | |
| □ NXT-35050 | 5" Polyurethane Light-Up | \$109 | | |
| □ NXT-34222 | 5" x 1.5" Poly, Plastic Mag | \$94 | | |
| □ NXT-42000 | 5" x 1.4" Poly, Aluminum Mag | \$99 | | |
| □ NXT-34666 | 6" x 1" Poly, Plastic Mag | No Charge | | |
| □ NXT-34999 | 6" x 1.5" Poly, Plastic Mag | \$94 | | |
| ■ NXT-44000 | 6" x 1.4" Poly, Aluminum Mag | \$99 | | |
| ■ NXT-34444 | 6" x 2" Poly with Tread | \$99 | | |
| ■ NXT-34667 | 6" x 1" Pneumatic Plastic Mag | \$83 | | |
| ■ NXT-34777 | 7" x 1" Poly, Plastic Mag | \$83 | | |
| □ NXT-34200 | 7" x 2" Polyurethane | \$127 | | |
| Note: All 7" casters not available with all 12" rear wheels. | | | | |
| | | | | |

REAR WHEELS, TIRES & ACCESSORIES

Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles. □ NXT-31524 Omit Rear Wheels and Axles

Note: All dual post armrests not available with 22" or 24" wheels.

| MAG WHEELS | Offili Rear Wheels and Axies | No Charge |
|----------------------------|--|-----------|
| □ NXT-31526 | 12" Low Profile Poly Mag | \$39 |
| | any foot/hand activated brake/wheel lock system. | 400 |
| Note: Not available with a | | |
| Note: All 12" wheels not a | available with all 7" casters. | |
| □ NXT-31523 | 12" Semi Pneumatic Plastic Mag | \$181 |
| Note: Not available with a | airless inserts. | |
| | available with all 7" casters. | |
| | 12" Pneumatic Plastic Mag | No Charge |
| | 19 12" Airless Inserts | \$83 |
| | available with all Handrims. | |
| | available with all 7" casters. | 0.4.4 |
| □ NXT-31530 | 16" Pneumatic Plastic Mag | \$44 |
| | 20 16" Airless Inserts I4 16" Aluminum Handrims | \$83 |
| | t available with plastic coated Handrims. | \$182 |
| □ NXT-31533 | 18" Pneumatic Plastic Mag | \$44 |
| | 21 18" Airless Inserts | \$83 |
| | 15 18" Aluminum Handrims | \$182 |
| | t available with plastic coated Handrims. | Ψ.02 |
| □ NXT-31527 | 20" Pneumatic Plastic Mag with Aluminum Handrims | No Charge |
| | 22 20" Airless Inserts | \$83 |
| O NXT-3101 | l6 20" Plastic ⊜ated Handrim Upgrade | \$128 |
| | | |

Select store location - AADL Numbers ☐ Calgary - AADL 00000 78000 003 ☐ Edmonton - AADL 00000 78000 001 ☐ Grande Prairie - AADL 0000 665538 001

No Charge

| FREEDOM NXI III ACC | cessories and Options | |
|---|---|-------------------|
| REAR WHEELS, TIRES & ACC | | |
| | 22" Pneumatic Plastic Mag with Aluminum Handrims 22" Airless Inserts | No Charge \$83 |
| | 22 Affiess fiserts 22" Plastic Coated Handrim Upgrade | \$128 |
| Note: All 22" & 24" wheels no | at available with short tilt bracket. | , |
| | ot available with dual post armrests. | No Chargo |
| | 24" Pneumatic Plastic Mag with Aluminum Handrims 24" Airless Inserts | No Charge \$83 |
| O NXT-31018 2 | 24" Plastic Coated Handrim Upgrade | \$128 |
| | t available with short tilt bracket. | |
| | t available with dual post armrests. Omit Handrims from 20", 22" or 24" Wheels | No Charge |
| SPOKE WHEELS | | |
| | 16" Grey Pneumatic Tire and Silver Spoke Aluminum Wheel | No Charge |
| Note: Rim, hub & spokes are | | φοο |
| | 16" Airless Inserts 16" Silver Aluminum Handrims | \$83 \$83 |
| | 16" Black Plastic Coated Handrims | \$83 |
| Note: 16" & 18" spoke wheels | | No Charge |
| □ NXT-31535 1 Note: Rim, hub & spokes are | 16" Black Pneumatic Tire and Black Spoke Aluminum Wheel | No Charge |
| | 16" Airless Inserts | \$83 |
| Note: 16" & 18" spoke wheels | | 000 |
| | 16" Black Aluminum Handrims 16" Black Plastic Coated Handrims | \$83 \$83 |
| Note: 16" & 18" spoke wheels | | Ç |
| EL NIVE 24527 | 40" Cook Brown of Time 9 Cilyan Cook a Alvaniana Miland | Na Obassa |
| □ NXT-31537 1 Note: Rim, hub & spokes are | 18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel | No Charge |
| · · · · · · · · · · · · · · · · · · · | 18" Airless Inserts | \$83 |
| | 18" Silver Aluminum Handrims | \$83 |
| O NXT-31036 1 Note: 16" & 18" spoke wheels | 18" Black Plastic Coated Handrims | \$83 |
| | 18" Black Pneumatic Tire & Black Spoke Aluminum Wheel | No Charge |
| Note: Rim, hub & spokes are | | 400 |
| | 18" Airless Inserts 18" Black Aluminum Handrims | \$83 \$83 |
| | 18" Black Plastic Coated Handrims | \$83 |
| Note: 16" & 18" spoke wheels | | |
| | 20" Pneumatic Spoke Wheel with Aluminum Handrims 20" Airless Inserts | \$220 \$83 |
| Note: Not available with foot/ | | φοσ |
| | heels not available plastic coated Handrims. | |
| | 22" Pneumatic Spoke Wheel with Aluminum Handrims 22" Airless Inserts | \$220 \$83 |
| Note: Not available with foot/l | | φοσ |
| | t available with short tilt bracket. | |
| | t available with dual post armrests. heels not available plastic coated Handrims. | |
| | 24" Pneumatic Spoke Wheel with Aluminum Handrims | \$220 |
| | 24" Airless Inserts | \$83 |
| Note: Not available with foot/l | | |
| | t available with short tilt bracket. It available with dual post armrests. | |
| Note: 20", 22" & 24" spoke wl | heels not available plastic coated Handrims. | |
| | Omit Handrims from 20", 22" or 24" Wheels | No Charge |
| REAR WHEELS, TIRES & ACC | Push to Lock Wheel Locks (Standard with NXT-30277) | No Charge |
| | Pull to Lock Wheel Locks | No Charge |
| | Wheel Lock Extension Handles | \$66 |
| | Foot Activated Wheel Lock (Drum) | No Charge |
| Note: Not available with 12" p | locks to be used with public transportation. poly plastic mag wheels. | |
| | Hand Activated Brake and Wheel Lock System (Drum) | No Charge |
| | locks to be used with public transportation. | |
| | nag wheels with low pro poly tires. n, speed, brake and lock wheels from strolling position. | |
| | Combination Dual Hand Activated Tilt, Brake and Wheel Lock System (Drum) | No Charge |
| | locks to be used with public transportation. | |
| Note: Not available with 12" p | noly plastic mag wheels. n, speed, brake and lock wheels and tilt frame from strolling position. | |
| REAR ANTI-TIPPERS | | |
| | Curb Assist Anti-Tip Assembly - Flip-Up, Adjustable | \$105 |
| | ip tube with curb assist feature. Standard Anti-Tip Assembly - Flip-Up, Adjustable | \$105 |
| SEAT PAN OPTIONS | Tanada a 7 ana 1 ap 7 tees and 3 y 1 ap 6 p , 7 tajastas a | ψ.00 |
| | Removable Solid Seat Pan - Bolteron | \$183 |
| | Removable Solid Seat Pan - Aluminum aises the effective seat to floor height by 5/8". | \$275 |
| | Solid Seat Pan adjusts 2" in width by mounting hardware to grow with NXT frame. | |
| Note: Removable solid seat p | pan includes 4-points of EZ-Lite mounting hardware. | |
| | " x 2" notches for growth adjustment. it frame dimensions unless specified below. | |
| | Seat Width" X Seat Depth" | |
| | at . an size. 4" extra depth maximum. | |
| | | |

□ Calgary - AADL 00000 78000 003
 □ Edmonton - AADL 00000 78000 001
 □ Grande Prairie - AADL 0000 665538 001

Select store location - AADL Numbers

| ACCESSORIES | | | | | |
|---|--|----------|--|--|--|
| □ NXT-23207 | Medical Necessities Pouch - Black Nylon | \$73 | | | |
| □ NXT-SFDI-9HP | Black Neoprene Frame Protectors - For Hangers | \$99 | | | |
| ☐ NXT-10911-2 | Black Neoprene Frame Protectors - For Back Posts | \$88 | | | |
| ☐ NXT-23210 | Flip-Down Rigidizing Bar | \$231 | | | |
| Note: Mounts to back post tubes to add rigidity to frame. | | | | | |
| Note: Does NOT adjust in | Note: Does NOT adjust in width to grow with frame. | | | | |
| ☐ NXT-10900 | Wheelchair Backpack | \$60 | | | |
| □ NXT-10905 | Wheelchair Canopy w/ Window and Storage Pocket - Flip-Back | \$138 | | | |
| Note: Canopy assembly includes removable mounting hardware. | | | | | |
| ANTERIOR SUPPORT SURFACES - TRAYS | | | | | |
| ☐ NXT-08-2000 | 3/8" Lexan Tray - Channel Bracket Hardware | \$292 | | | |
| ☐ NXT-08-6000 | 3/8" Lexan Tray - Adjustable Knob Hardware & Ribs | \$236 | | | |
| □ NXT-08-13000 | 1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs | \$197 | | | |
| ☐ NXT-08-16000 | 1/4" Lexan Tray - Hook & Loop Straps | \$114 | | | |
| ■ NXT-10917 | Wheelchair Tool Kit | Standard | | | |
| SEATING | | | | | |

SEA

Note:

Note:

Note:

| | □ NXT-08-13000 □ NXT-08-16000 | 1/4" Lexan Tray - Hook & Loop Straps Wheelshair Tray - Hook & Loop Straps | \$197 \$114 | | |
|-------------------------------------|---|--|----------------|--|--|
| | ■ NXT-10917 | Wheelchair Tool Kit | Standard | | |
| ATII | | ly of the chair | | | |
| | seating will ship separate t all seating options are lis | sted. See the full Matrx seating order form for more options. | | | |
| | | rd sizes only unless otherwise noted. | | | |
| | BACK SUPPORTS | | | | |
| MATRX MINI - PEDIATRIC POSTURE BACK | | | | | |
| | ☐ MINIXXXX-S | Black Startex Fabric Cover | \$845 | | |
| | □ MINIXXXX-K | Kid*a*bra Meshtex Fabric Cover | \$845 | | |
| | ☐ MINIXXXX-B | Blue Blizzard Startex Fabric Cover | \$845 \$845 | | |
| | ☐ MINIXXXX-P | Pink Lollipop Startex Fabric Cover | \$845 | | |
| | | ECLINE BACK WITH SUPERSOFT FOAM | | | |
| | | layer of high resiliency foam for added comfort and immersion. | # 000 | | |
| | □ PBEXXXX-TR | Elite TR Back Wingship Slippover Add On for Elite TR Back (Stock sizes only) | \$866 \$150 | | |
| | □ PBEWSCXXXX-IR | Wipeable Slipcover Add-On for Elite TR Back (Stock sizes only) | \$150 | | |
| | ELITE - POSITIONIN | | | | |
| | □ PBEXXXX | PB Elite Back | \$766 | | |
| | ☐ PBEWSCXXXX | Wipeable Slipcover Add-On for Elite Back (Stock sizes only) | \$150 | | |
| | | | | | |
| | CUSHIONS | | | | |
| | KID*A*BRA VI - PED ☐ KAXXXX-VI | IATRIC SKIN PROTECTION AND POSITIONING FOAM CUSHION Kid*a*bra Meshtex Fabric Cover | \$545 | | |
| | | C POSTURE SEAT FOAM CUSHION | | | |
| | □ PSXXXX | Reversible Black Startex Cover | \$569 | | |
| | □ PSXXXX-K | Kid*a*bra Startex Fabric Cover (Non-reversible) | \$569 | | |
| | □ PSXXXX-B | Blue Blizzard Startex Fabric Cover (Non-reversible) | \$569 | | |
| | □ PS XXX-P | Pink Lollipop Startex Fabric Cover (Non-reversible) | \$569 | | |
| | LIBRA - SKIN PROT | ECTION AND POSITIONING FOAM AND FLUID CUSHION | | | |
| | □ LCXXXX | Startex Cover, Fabric Side Out <= 20W x 20D | \$745 | | |
| | □ LCXXXX-IC | Infection Control Cover <= to 20W x 20D | \$798 | | |
| | | | | | |
| | ACCESSORIES | | | | |
| | BASIC HEADREST | | | | |
| | □ 5500071 | Comfort Headrest Assembly w/ Multi Axis Hardware | \$375 | | |
| | ELAN HEADREST P | ADS - HARDWARE SOLD SEPARATELY | | | |
| | □ ESP6 | Standard Pad - Startex - tailored, removable cover. 6" x 3.5" | \$310 | | |
| | □ ESP10 | Standard Pad - Startex - tailored, removable cover. 10" x 5" | \$310 | | |
| | □ ESP14 | Standard Pad - Startex - tailored, removable cover. 14" x 5" | \$310 | | |
| | ☐ ESP6-IC | Standard Pad - Infection Control - tailored, removable cover. 6" x 3.5" | \$335 | | |
| | ☐ ESP10-IC | Standard Pad - Infection Control - tailored, removable cover. 10" x 5" | \$335 | | |
| | ☐ ESP14-IC | Standard Pad - Infection Control - tailored, removable cover. 14" x 5" | \$335 | | |
| | □ EOP9 | Occipital Pad - Startex - tailored, removable cover. 9" x 7" | \$341 | | |
| | □ EOP12 | Occipital Pad - Startex - tailored, removable cover. 12" x 8" | \$341 | | |
| | ☐ EOP9-IC | Occipital Pad - Infection Control - tailored, removable cover. 9" x 7" | \$366 | | |
| | ☐ EOP12-IC | Occipital Pad - Infection Control - tailored, removable cover. 12" x 8" | \$366 | | |
| | ☐ E4POINT | 4 Point Postioning Pad - Startex - tailored, removable cover | \$448 | | |
| | ☐ E4POINT-IC | 4 Point Postioning Pad - Infection Control - tailored, removable cover | \$473 | | |
| | ELAN HEADREST H | ARDWARE - PADS SOLD SEPARATELY | | | |
| | | Standard Flan Hardware 15" Vertical Post | \$266 | | |

Standard Elan Hardware - 15" Vertical Post Mini Elan Hardware - 11" Vertical Post



Kid*a*bra Meshtex Fabric



Blue Blizzard Startex Fabric



Pink Lollipop Startex Fabric

NOTES:



□ MEHW

□ MEMH

Select store location - AADL Numbers

- ☐ Calgary AADL 00000 78000 003
- ☐ Edmonton AADL 00000 78000 001
- Grande Prairie AADL 0000 665538 001

\$266

\$254