Cruiser X4 Wheelchair

HCPCS: K0004-1,K0004



• Weighs 34 pounds or less when properly equipped in the basic configuration (Excluding front rigging)

• Easy to propel and transport

drive

• Frame style eliminates seat guides and allows for custom backs and accessories

• Flip-back, removable arm is height adjustable from top of seat rail 10" to 14" and features a front mounted easy-to-use lever release mechanism (Figure A - C)

· Padded armrests provide added patient comfort

🗇 DeVilbiss

• Height adjustable segmented back

• Aluminum and Carbon steel blended frame with combination powder coated and silver vein finish is attractive, chip resistant and easy-to-maintain

- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Patented adjustable depth seat 16" and 18" on chair depths (Figure D F)
- 3 position quick release rear axles mounted on Composite Mag-style wheels with composite hand rims are lightweight and maintenance free (Figure G)
- Polyurethane rear tires and composite casters offer superior performance, smoother ride and are lightweight
- Dual over center cross links provide added rigidity (Figure H)
- · Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- · Composite footplates with heel loops are durable and lightweight
- Precision sealed wheel bearings through-out ensure long lasting and reliability
- 8" front casters have 3 height adjustments and angle adjustments (Figure I)
- Triple rear axle positions provides easy transition of seat height from standard to hemi level
- Comes with push-to-lock wheel locks

PRODUCTS	Description
CX416ADDA-ELR	CruiserX4, 16" Adj DeskArm, ElvLg, 1/ea
CX416ADDA-SF	CruiserX4, 16" Adj DeskArm, SwgFt, 1/ea
CX416ADFA-ELR	CruiserX4, 16" Adj FullArm, ElvLg, 1/ea
CX416ADFA-SF	CruiserX4, 16" Adj FullArm, SwgFt, 1/ea
CX418ADDA-ELR	CruiserX4, 18" Adj DeskArm, ElvLg, 1/ea
CX418ADDA-SF	CruiserX4, 18" Adj DeskArm, SwgFt, 1/ea
CX418ADFA-ELR	CruiserX4, 18" Adj FullArm, ElvLg, 1/ea
CX418ADFA-SF	CruiserX4, 18" Adj FullArm, SwgFt, 1/ea
CX420ADDA-ELR	CruiserX4, 20" Adj DeskArm, ElvLg, 1/ea
CX420ADDA-SF	CruiserX4, 20" Adj DeskArm, SwgFt, 1/ea
CX420ADFA-ELR	CruiserX4, 20" Adj FullArm, ElvLg, 1/ea
CX420ADFA-SF	CruiserX4, 20" Adj FullArm, SwgFt, 1/ea

SPECIFICATIONS	Description
Back Height	16" - 17.5"
Carton Shipping Weight	46 lbs 48 lbs 50 lbs
Height	36"
Length with Riggings	43"
Product Weight without Riggings (16, 18, 20	32 lbs 34 lbs 36 lbs
)	
Seat Depth	16"/18" On Chair
Seat To Floor Height	17" 18" 19"
Seat Width	16" 18" 20"
Warranty	Limited Lifetime
Warranty on Frame	Limited Lifetime
Weight Capacity	300 lbs
Width Closed	12"
Width Open	23" 25" 27"
ACCESSORIES	Description
STDSSHK-K3	Super Hemi Kit, Cruiser 3 & X4, 1/ea, EA
STDS851	Seat Belt, Velcro Close, 1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS850	Seat Belt, Clasp Type,1/ea, EA
STD5834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
WASL	Limb Support,Left,1/cs, CS
WASE	Limb Support, Right, 1/cs, CS
STDS601	Brake Extension,6",1/ea, EA
STD5001 STD51002	Caster,8in Sentra,Viper,Astaire,1/ea, EA
STDS1002 STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
5125000	
PARTS	Description
STDS4SC4X-03	LeftFrntCasterHousingW/Bearings&HrdwrCX4, EA
STDS4SC4X-04	RightFrntCasterHousingW/Bearings&HrdwCX4, EA
STDS4SC4X-17	Fork, Cruiser 4X, 1/ea, EA
STDS3J24AFL	Adj. Full Arm, Cruiser 3 / X4, Left, 1e, EA
STDS3J24AFR	Adj. Full Arm, Cruiser 3 / X4,Right,1ea, EA
STDS3J24AFR-UFR	UpprSectOnlyW/Pad,AdjFullArm,Rt,Crsr3/4, EA
STDS3J24AL-UDL	UpprSectOnlyW/Pad,AdjDskArm,Lft,Crsr3/4, EA
STDS4SC4X-08	HeightAdjustmentLever,Left,CruiserX4,1e, EA
STDS4SC4X-09	HeightAdjustmentLever,Rght,CruiserX4,1e, EA
STDS4SC4X-15	PushHandleW/Grip&AttachHrdwr,CruiserX4, EA
STDS4SC4X-10	Back Upholstery, 16", CruiserX4, 1/ea, EA
STDS4SC4X-11	Back Upholstery, 18", CruiserX4, 1/ea, EA
STDS4SC4X-12	Back Upholstery, 20", CruiserX4, 1/ea, EA
STDS4SC4X-01	RearWheelComplW/Bearings&Hndrm,Cruiser4, EA
STDS4SC4X-02	Quick Release Rear Axle, CruiserX4, 1/e, EA
STDS4SC4X-16	RearWheelAxleReceptacle, CruiserX4, 1/e, EA
STDS4SC4X-05	Side Frame, Left, CruiserX4, 1/ea, EA
STDS4SC4X-06	Side Frame, Right, CruiserX4, 1/ea, EA
STDS4S1069	HangerPins,CS/SilverSport/Cruiser, 4/pkg, PKG
STDS3J24X6-AD	Cross Frame,Cruiser3,Adj Depth5D,16",1ea, EA