# PRODUCT DATA SHEET



# THE WORLD'S NO.1 STAIRWAY EVACUATION CHAIR

# Evac+Chair® Power 800

THIS ELECTRIC/BATTERY POWERED STAIR CLIMBER ALLOWS FOR EASY AND SAFE STAIR ASCENT & DESCENT

Can be used in emergency evacuations, power outages, and daily patient transport to safely and easily move a passenger up and down a flight of stairs.

#### **Features**

#### 1 Height adjustable head cushion

Headrest with padded cushion for maximum comfort for the passenger during transit

## 2 Easy to use control panel

Single speed control panel with on/off power ar 1 ascent and descent buttons

# 3 Ergonomic adjustable handle

For ease of deployment and compact storage

### 4 Additional operator handles

For added stability when using additional operators during transportation

# 5 Extra long quick release four point harness

For increased passenger safety

# 6 Rechargeable lithium Ion battery pack

Fully charged in 1.5 hours (Equivalent to Going Up and Down a 5 Story building flight of stairs 30 times)

#### 7 High strength aluminium frame

Increased strength and robustness for continued use

## 8 LED safety light and photoluminescent branding

For ease of use in situations of poor visibility

#### 9 Adjustable stair track


Tension control system for perfect transportation with tool included

# 10 Lockable rear wheels


Can be applied in situations where the chair needs to remain in a stationary position

Dimensions	Height	Width	Depth	Weight	Weight Capacity
Folded	45"	19.5"	10.6"	66lbs	350lbs

9









# **Includes**

- Battery Charger (US/EU/UK)
- Rechargeable Battery
- 2 Keys
- User Guide
- Passenger Headrest
- Passenger Safety Straps
- Photoluminecent Wall Sign
- Wall Hooks
- **Dust Cover**
- Manufacturers Warranty

#### **Benefits**

- FDA Registered
- CE & ISO 9001 Certificated
- Robust and compact in size
- Provides a quick, easy and safe emergency evacuation
- Easy to use
- One person operation
- No heavy lifting required
- Enables small attendants to easily move passengers three times their size
- Easy to open, ready to use in seconds
- Folds for compact storage
- Flame retardant materials
- Wipeable materials and anti-microbial paint for infection control
- Exceeds fire safety requirements and emergency response guidelines established to support the Americans with Disabilities Act

# VIDEO LINK: evac-chair.ca/watch-video/

For over 35 years EVAC+CHAIR® is America's Original and World's Number One Selling Evacuation Chair. Don't Risk Lives, Beware of Poor Imitations and Inferior Substitutions. There is only one EVAC+CHAIR®. Unsurpassed Quality and Workmanship, Successful in Saving Lives in 9/11.

For further information please contact:

