

2 in 1 TRANSPORT CHAIR Series

THE TRANSPORT CHAIR COMBINES BOTH WALKER AND WHEELCHAIR INTO ONE EASY-TO-USE UNIT.

The center-folding design and folding footrest allows for easy storage. Back rest swings back for use as a wheelchair and forward for use as a walker.

FEATURES

BACK SUPPORT Extra wide back strap for increased support

TOTE BAG Zippered tote bag for privacy

CENTRE-FOLDING Small footprint when not in use

WHEELS 8"x1.5" tires with dual bearings

COLORS Warm Grey



SPECIFICATIONS

20 lbs
33.5" - 38.5" Handle Height
24" Width x 39" Length (Unfolded)
13.5" Width x 28.5" Length (Folded)
300 lbs Weight Capacity



*Net Weight of walker without accessories, backstrap and folding basket