

# PIPER Series



## SPECIFICATIONS

### SUPER MINI

12 lbs\*  
16.5" Seat Width  
17" Seat Height  
21" Width x 26" Length  
29" - 33" Handle Height  
275 lbs Weight Capacity



### MINI

13 lbs\*  
18" Seat Width  
19" Seat Height  
21" Width x 26" Length  
30" - 34" Handle Height  
275 lbs Weight Capacity



### REGULAR

13 lbs\*  
18" Seat Width  
21" Seat Height  
21" Width x 26" Length  
31" - 36" Handle Height  
275 lbs Weight Capacity



### TALL

13 lbs\*  
18" Seat Width  
23" Seat Height  
23" Width x 26" Length  
35" - 39" Handle Height  
275 lbs Weight Capacity

THE BEST SELLING WALKER IN ITS CLASS, WAS DESIGNED FOR THOSE WHO NEED A LIGHTWEIGHT YET SOLID AND STABLE PLATFORM.

It's the one of the strongest aluminum walkers on the market today and is available in four sizes. Piper Super Mini with a 17" seat height weighing just 12 lbs to the Piper Tall with a 23" seat height weighing just 13 lbs.

## FEATURES

**ALUMINUM FRAME** Superb stability, and style

**SEAT** Solid plastic base with removable fitted waterproof cushion

**MESH BASKET** Removable

**WHEELS** 8" 40 mm wide, Solid grip / Longer lasting polyurethane, Suitable for indoor & outdoor use

**COLORS** Vibrant Red, Black



## OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Deluxe Back Strap Pad
- Oxygen Tank Holder
- Mono Brakes
- Tote Bag (fits inside basket)
- Tray - Seat Mounted

\*Net Weight of walker without accessories, backstrap and folding basket