



Folding Junior Walker, Two Button with 5" Wheels

HCPCS: E0143-1,E0143

- Easy push-button release mechanisms
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter construction ensures maximum strength
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces
- Vinyl-contoured hand grip
- Comes with wheels
- Limited Lifetime Warranty

| PRODUCTS | Description |
|------------------------|---|
| 10247-4 | Walker,Steel, Junior 2Button 5" Whs,4/CS |
| SPECIFICATIONS | Description |
| Carton Shipping Weight | 30.64 lbs |
| Depth (opened at base) | 18" |
| Height | 25" - 32" |
| Overall Width | 24" |
| Warranty | Limited Lifetime |
| Weight | 8 lbs |
| Weight Capacity | 350 lbs |
| Width | 17" inside hand grip 20" inside back legs |