

Maple Leaf WHEELCHAIR



The only limitation is your imagination

JR TILT Commode



PEDIATRIC TILT
COMMODOE



www.mapleleafwheelchair.ca

JR **TILT**Commode



◆ **GROWING FRAME**

The JR TILT COMMODOE is Maple Leaf Wheelchair's pediatric tilting commode that can be customized to a wide range of specifications, fits over standard toilets, and comes standard with a growing frame increasing in size as the child does!

STANDARD FEATURES

FRAME & UPHOLSTERY

- ◆ 14" - 16" growing frame width
- ◆ 14" - 16" growing frame depth
- ◆ 21" or 23" seat to floor height
- ◆ 16" back height
- ◆ Soft seat with plastic bucket
- ◆ 4-way adjustable contoured multi axis headrest
- ◆ Nickel plated epoxy coated frame
- ◆ Strong and flexible adjustable water proof mesh back
- ◆ Laterals

TILT

- ◆ Stainless steel inner cable with PTFE lined conduit
- ◆ Gas struts for tilting
- ◆ Anti-tippers

ARMRESTS & FOOTRESTS

- ◆ Durable urethane arm pads
- ◆ Flip back armrests
- ◆ Swing away footrests with composite footplates

WHEELS

- ◆ 5" rust resistant nylon casters
 - ◆ (3) - total lock casters
 - ◆ (1) - directional lock caster

ACCESSORIES & OPTIONS

- ◆ Custom seat width
- ◆ Custom seat depth
- ◆ Custom seat heights
- ◆ Custom back heights
- ◆ Custom seating
- ◆ Growing kits
- ◆ 24" rear wheels
- ◆ Extra side laterals

*see full list of accessories on the JR TILT COMMODOE order form

MODEL SPECIFICATIONS

- ◆ 5 year warranty on frame
- ◆ 2 year warranty on additional components

Dealer Information

6540 Tomken Road
Mississauga ON L5T 2E9
T 905.564.2250 F 905.564.5424
1.800.362.0873

www.mapleleafwheelchair.ca

