

Permobil

Heavy duty products



permobil



Up to 8" of seat elevation

Base width: 26"
Base length: 38"

Permobil M300 Corpus® HD

The M300 Corpus HD takes all the things you love about Permobil's mid-wheel drive and amplifies them. Outfitted with stronger actuators, more powerful motors, and a reinforced Corpus seating system, it's built to withstand the rigors of an active heavy duty lifestyle. Designed specifically for users up to 450 lb, the suspension has been strengthened to provide you with more stability and completed with elegant features like fully adjustable legrests and up to four power seat functions.



Weight capacity— 450 lb


Specifications

Max user weight	450 lb
Max speed	5 mph
Range	13 mi (Grp 34) 20 mi (Grp 24)
Base width	26"
Base length	38"
Seat-to-floor height	18", 19", 20" (elevator) 18", 19", 20" (fixed)
Base turning radius	22"
Ground clearance	3"
Recommended battery type	Group 34 or 24
Weight including batteries	425 lb (Group 34) 447 lb (Group 24)
Drive electronics	R-Net 120A

Optional seat elevator	8"
Tilt options	Posterior: 0° – 45°
Backrest angle - power recline	85° – 150°
Legrest angle - power	85° – 175°
Seat depth	19" – 25" (by 1" incr.)
Seat width	19", 21", 23", 25"
Backrest height	21", 24" – 29" (by 1" incr.)

* Actual driving range will vary based on driving and battery conditions.
All specifications may vary based upon configuration.
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.

Color options



Permobil M3 Corpus

Meticulously engineered to help you maintain better positioning, feel more secure and reduce driving fatigue so you can conquer daily activities with ease. Its compact base offers maneuverability in tight spaces without sacrificing performance. Featuring the innovative FlexLink suspension geometry with independent caster swing-arms helping to keep you stable and in control as you traverse the obstacles in your world at up to 6 mph.




Specifications

Max user weight	330 lb
Max speed	6 mph
Range	16 mi (Grp 34) 20 mi (Grp 24)
Base width	24"
Base length	35.5"
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	20.5"
Ground clearance	3"
Recommended battery type	Group 34 or 24
Weight including batteries	385 lb (Group 34) 405 lb (Group 24)
Drive electronics	R-Net 120A
Optional seat elevator	12"

Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° (up to 300 lb)
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

** Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

Color options



Chili Sunset Canary Spring Aqua Sky Violet Flamingo Arctic Midnight

Premium options



Gold Platinum Bronze

Permobil M5 Corpus

Take control with high-performance motors enabling a top speed of 7.5 mph and exceptional low-end torque that lets you pivot like a pro. Experience unsurpassed stability and ride comfort with FlexLink Pro independent drive wheel suspension. The M5 Corpus will take you to new heights with increased stability and full suspension in seat elevation.

Up to 14" of ActiveHeight™ & 20° of ActiveReach™

Base width: 25"
Base length: 37.5"



Weight capacity— 330 lb

Specifications

Max user weight	330 lb
Max speed	7.5 mph
Range	20 mi
Base width	25"
Base length	37.5"
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	21.5"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	423 lb
Drive electronics	R-Net 120A
Optional seat elevator	14"


Tilt options	Posterior: 0° – 50°
	ActiveReach: 5°, 10° & 20° (up to 300 lb)
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

* Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.

Color options



Premium options



Permobil F3 Corpus

The F3 Corpus with a 6 mph top-end speed features optimized battery weight distribution, lower center of mass, and fully adjustable ComfortRide suspension that improve stability while traversing inclines or declines. Experience the enhanced functionality that combines outstanding chassis performance with Permobil's industry-leading Corpus seating system.



Specifications

Max user weight	330 lb	Drive electronics	R-Net 120A
Max speed	6 mph 5 mph <i>without support wheels</i>	Optional seat elevator	12"
Range	16 mi (Grp 34) 20 mi (Grp 24)	Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20° & 30°
Base width	24"	Backrest angle - power recline	85° – 180°
Base length	40" with support wheels 36" without support wheels	Legrest angle - power	85° – 180°
Seat-to-floor height	17.5", 18.5", 19.5"	Seat depth	14" – 22" (1" incr.)
Base turning radius	26.5"	Seat width	17" – 23" (2" incr.)
Ground clearance	3"	Backrest height	20", 23" – 28" (1" incr.)
Recommended battery type	Group 34 or 24	<i>* Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.</i>	
Weight including batteries	386 lb (Group 34) 407 lb (Group 24)		


Color options



Premium options

Permobil F5 Corpus

Power meets performance on the F5 Corpus with a 7.5 mph top-end speed and upgraded ComfortRide Pro suspension. All-wheel independent, fully adjustable, oil dampened shocks provide the ultimate in ride comfort and reduced vibration while driving in all seat positions, including an unprecedented 14" ActiveHeight seat elevation. Climb and reach higher with increased obstacle negotiation, exceptional motor torque, and enhanced stability allowing up to 45° of ActiveReach forward tilt functionality.




Specifications

Max user weight	330 lb
Max speed	7.5 mph 6 mph without support wheels
Range	20 mi
Base width	25.5"
Base length	42" with support wheels 38" without support wheels
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	30"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	412 lb
Drive electronics	R-Net 120A

Optional seat elevator	14"
Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20°, & 45°
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)


** Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

Color options



Chili Sunset Canary Spring Aqua Sky Violet Flamingo Arctic Midnight

Premium options



Gold Platinum Bronze

TiLite Aero T

Our most advanced aluminum chair is uniquely modeled and handmade to each user's specifications. The dual-tube frame of the Aero T provides superior rigidity, while minimizing weight and maximizing performance. Thanks to the TiFit process, you get an optimized fit plus adjustability that meets your needs.



Note: The ROHO® AGILITY™ Active Carbon back supports up to 265 lb.




Specifications

Max user weight	265 lb standard 350 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium dual-tube rigid, TiFit, full adjustability	Center of gravity adjustment	5"
Total chair width*	23"	Front seat width (seat taper)	0", 1", 2", 3"
Seat width	12" – 20"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	12" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	15" – 21"	Transit options	Q-Straint lockdown approved
Rear seat height	12.5" – 21"	Transport weight*	12 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


Frame colors

High gloss paint




TiLite red Neon pink High gloss black

Candy paint




Ruby red Azure blue Electric plum

Metallic paint



Tangerine metallic Acid green metallic Ocean blue metallic

Matte textured paint




Flame orange OD green Hyper blue

Matte paint



Titanium grey Black lava Midnight blue pearl


Pearlescent paint



White pearl

Component colors


Anodized



Red Orange Gold



Green Blue Purple



Black



TiLite TRA

With the power and style of its titanium dual-tube frame, plus TiFit precision for an exact fit, the TRA provides performance and adjustability made just for you.



Weight capacity— 350 lb


Specifications

Max user weight	265 lb standard 350 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium dual-tube rigid, TiFit, full adjustability	Center of gravity adjustment	5"
Total chair width*	23"	Front seat width (seat taper)	0", 1", 2", 3"
Seat width	12" – 20"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	12" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	15" – 21"	Transit options	Q-Strait lockdown approved
Rear seat height	12.5" – 21"	Transport weight*	11.1 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


Frame colors

High gloss paint




TiLite red Neon pink High gloss black

Candy paint




Ruby red Azure blue Electric plum

Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

Matte textured paint




Flame orange OD green Hyper blue

Matte paint




Titanium grey Black lava Midnight blue pearl

Pearlescent paint



White pearl Satin Polished


Natural titanium



Green Blue Purple

Component colors

Anodized



Red Orange Gold
Green Blue Purple
Black

TiLite ZRA

With full adjustability and a minimalist mono-tube frame, the ZRA dishes out style without boundaries. With its TiFit titanium frame, you have the complete package: strength, style, and the ability to set up and configure your chair just for you.



Weight capacity— 300 lb




Specifications

Max user weight	265 lb standard 300 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium mono-tube rigid, TiFit, full adjustability	Center of gravity adjustment	6"
Total chair width*	23"	Front seat width (seat taper)	0", 1.5", 2.5"
Seat width	10" – 20"	Precision Lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	10" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	16" – 21"	Transit options	Transit tie-down
Rear seat height	12.5" – 20"	Transport weight*	10.7 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


Frame colors

High gloss paint




TiLite red Neon pink High gloss black

Candy paint




Ruby red Azure blue Electric plum

Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

Matte textured paint




Flame orange OD green Hyper blue

Matte paint



Titanium grey Black lava Midnight blue pearl

Pearlescent paint




White pearl Satin Polished

Natural titanium

Component colors

Anodized



Red Orange Gold
Green Blue Purple
Black



TiLite ZR

With the ZR's minimalist mono-tube frame, less truly is more. The clean aesthetic of the ZR turns heads while the titanium frame makes sure you have the strength and vibration-dampening ride to go with your beauty. Style and strength to take on whatever life throws at you.



Weight capacity— 300 lb


Specifications

Max user weight	265 lb standard 350 lb heavy duty	Center of gravity adjustment	2.25"
Frame features	Titanium dual-tube rigid, TiFit, limited adjustability	Front seat width (seat taper)	0", 1.5", 2.5"
Total chair width*	23"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat width	10" – 20"	Fixed seat back	Fixed height: 8" – 20" Fixed angle: 80° – 101°
Seat depth	10" – 20"	Ergonomic seat	4" – 10"
Front seat height	16" – 21"	Transit options	Transport bracket
Rear seat height	12.5" – 20"	Transport weight*	9 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components
Camber	0°, 2°, 4°, 6°, 8°, or 12°		

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


Frame colors

High gloss paint




TiLite red Neon pink High gloss black

Candy paint




Ruby red Azure blue Electric plum

Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

Matte textured paint




Flame orange OD green Hyper blue

Matte paint




Titanium grey Black lava Midnight blue pearl

Pearlescent paint



White pearl Satin Polished


Natural titanium



Black

Component colors

Anodized



Red Orange Gold
Green Blue Purple
Black

TiLite Aero X

The Aero X features an innovative modular design with clean modern lines and configurability. Featuring multiple configurations that include heavy duty, hemi-height, and amputee options. Add superior rolling dynamics and full adjustability, and you get a chair that excels in both performance and price.




Specifications

Max user weight	265 lb standard 350 lb heavy duty	Front frame angle	Fixed front: 70°, 80°, or 85° Swing-away: 70°, 80°, or 85°
Frame features	Aluminum box folding frame, modular design, full adjustability	Camber	0°, 2°, 4°, or 6°
Total chair width folded*	10.5"	Center of gravity adjustment	6"
Total chair width*	22.5"	Fixed seat back	Fixed height: 11" – 20"
Seat width	12" – 22"	Transit options	Transit tie-down, transport bracket
Seat depth	12" – 20"	Transport weight*	Fixed front: 15.6 lb Swing-away: 18 lb
Front seat height	13.5" – 21.5"	Warranty	Five years on frame One year on TiLite components
Rear seat height	13" – 20"		

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


Frame colors

High gloss paint




TiLite red Neon pink High gloss black

Candy paint




Ruby red Azure blue Electric plum

Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

Matte textured paint




Flame orange OD green Hyper blue

Matte paint



Titanium grey Black lava Midnight blue pearl

Pearlescent paint



White pearl



TiLite 2GX

The TiLite 2GX combines the advantages of a modular folding frame with the high strength and vibration-dampening characteristics of titanium. Add in the Tru-Fit full adjustability that meets your needs and you are ready to roll in a chair designed to empower you.




Specifications

Max user weight	265 lb standard 350 lb heavy duty	Front frame angle	Fixed front: 70°, 80°, or 85° Swing-away: 70°, 80°, or 85°
Frame features	Titanium box folding frame, modular design, full adjustability	Camber	0°, 2°, 4°, or 6°
Total chair width folded*	10.75"	Center of gravity adjustment	6"
Total chair width*	23"	Fixed seat back	Fixed height: 11" – 20"
Seat width	18" – 22"	Transit options	Transport bracket
Seat depth	16" – 20"	Transport weight*	Fixed front: 13.1 lb Swing-away: 15.5 lb
Front seat height	13.5" – 21"	Warranty	Five years on frame One year on TiLite components
Rear seat height	13" – 20"		

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


Frame colors

High gloss paint




TiLite red Neon pink High gloss black

Candy paint




Ruby red Azure blue Electric plum

Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

Matte textured paint




Flame orange OD green Hyper blue

Matte paint




Titanium grey Black lava Midnight blue pearl

Pearlescent paint



White pearl

Natural titanium



Satin Polished



SmartDrive MX2+

SmartDrive's intuitive controls add a power assist for a variety of manual wheelchair scenarios like cruising up hills or over thick carpets. Its lightweight design attaches quickly and its unique OmniWheel design moves with you as you pivot, turn, wheelie, and navigate curbs in style. SmartDrive gives you freedom to do more.



Weight capacity— 331 lb

Specifications

Max user weight	30 – 331 lb
Product dimensions	Length: 15.3" Width: 5.7" Height: 9.5"
Product weight	13.2 lb drive unit
Driving range*	12.3 mi level/no resistance
Motor power	250 W brushless DC
Maximum speed	5.5 mph
Operating temp	-13°F – 122°F
Chair type compatibility	Rigid, folding, tilt-in-space, one-arm drive, and standing

Drive wheel diameter	22", 24", 25", 26"
Warranty	2-Year limited warranty
Bluetooth® LE	FCC ID QOQBT113 IC ID 5123A-BGTBLE113

**Performance measurements will vary based on user and chair characteristics, driving and battery conditions.*

The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.

A wheelchair's specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.

Cushions & back supports



ROHO Contour Select

No weight limit when properly sized to cushion and individual.



ROHO® Hybrid Elite® SR® Single Compartment cushion*

Ideal for individuals who May require accommodation of pelvic asymmetry, up to approx. 2"
May require pelvic and lower extremity positioning and a more stable environment for transfers

Dimensions	Width: 13.75" – 23.75"
	Depth: 14.75" – 22.75"
	Height: 4.25" without load
Weight limit	Available up to 500 lb
Medicare codes	E2622, E2623



ROHO Hybrid Elite Dual Compartment cushion*

May require correction of pelvic asymmetry (medial lateral), up to approx. 2"
May require a more stable environment for transfers

Dimensions	Width: 13.75" – 23.75"
	Depth: 14.75" – 22.75"
	Height: 4.25" without load
Weight limit	Available up to 500 lb
Medicare codes	E2622, E2623



ROHO CONTOUR SELECT® cushion

May require accommodation of pelvic asymmetry, up to approx. 1"
May need correction of pelvic asymmetry (medial lateral AND anterior posterior)
May require different positioning of the left and right femur

Dimensions	Width: 14.5" – 24.25"
	Depth: 15.25" – 22"
	Height: 2.25" – 4.25"
Weight limit	No weight limit on properly sized cushions
Medicare codes	E2624, E2625

*Hybrid Elite images do not show the cushion cover, which is required for use.



ROHO QUADTRO SELECT® HIGH PROFILE® cushion

Ideal for individuals who May be at risk for or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
May require accommodation or correction of medial lateral and/or anterior posterior pelvic asymmetries, up to approx. 3"
May require different positioning of the left and right femur

Dimensions	Width: 11" – 26.25"
	Depth: 11.75" – 23.75"
	Height: 4.25"
Weight limit	No weight limit on properly sized cushions
Medicare codes	E2624, E2625



ROHO QUADTRO SELECT® MID PROFILE® cushion

May be at risk for or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
May require accommodation or correction of medial lateral and/or anterior posterior pelvic asymmetries, up to approx. 2"
May require different positioning of the left and right femur

Dimensions	Width: 11" – 26.25"
	Depth: 11.75" – 23.75"
	Height: 3.25"
Weight limit	No weight limit on properly sized cushions
Medicare codes	E2624, E2625



ROHO ENHANCER® cushion

May require accommodation of pelvic asymmetry, up to approx. 1"
Require positioning of the pelvis and lower extremities
May need increased lateral support
May need to partially offload the ischial tuberosities

Dimensions	Width: 14.75" – 20"
	Depth: 14.75" – 20"
	Height: 2.25" – 4.25"
Weight limit	No weight limit on properly sized cushions
Medicare codes	E2624, E2625



ROHO HIGH PROFILE Dual Compartment cushion

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
 May require accommodation of pelvic asymmetry, up to approx. 3"
 May require correction of pelvic asymmetry (medial lateral and/or anterior posterior)



ROHO HIGH PROFILE Single Compartment cushion

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
 May require accommodation of pelvic asymmetry, up to approx. 3"
 May require maximum amount of immersion and envelopment



ROHO MID PROFILE Single Compartment cushion

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
 May require accommodation of pelvic asymmetry, up to approx. 2"
 May require immersion and envelopment

Ideal for individuals who

Dimensions

Weight limit

Medicare codes

Width: 14.75" – 20"	Width: 11.5" – 26.75"	Width: 11.5" – 26.75"
Depth: 14.75" – 21.75"	Depth: 11.5" – 21.75"	Depth: 11.5" – 21.75"
Height: 4.25"	Height: 4.25"	Height: 3.25"
No weight limit on properly sized cushions	No weight limit on properly sized cushions	No weight limit on properly sized cushions
E2624, E2625	E2622, E2623	E2622, E2623



Comfort Support Pro 7 Series cushion

May need relief in the coccyx area
 Requires/wants cushion to strap to chair



Comfort Saddle 7 Series cushion

May need postural alignment
 May need relief in the coccyx area
 May require minimal to moderate skin protection
 Requires/wants cushion to strap to chairs

Ideal for individuals who

Dimensions

Weight limit

Medicare codes

Width: 22" – 30"	Width: 22" – 30"
Depth: 18" – 20"	Depth: 18" – 20"
750 lb	750 lb
E2602	E2608



Comfort Embrace cushion

May need postural alignment
 May need support to the greater trochanters



Comfort Embrace cushion with ATI

May need postural alignment
 May need support to the greater trochanters
 May need dynamic positioning



Comfort Ascent cushion

May need support for the greater trochanters
 May require mild-moderate skin protection

Ideal for individuals who

Dimensions

Weight limit

Medicare codes

Width: 18" – 24"	Width: 18" – 24"	Width: 16" – 24"
Depth: 14" – 24"	Depth: 14" – 24"	Depth: 16" – 20"
Height: 4" for Zero elevation, 3.75" for Anti-Thrust	Height: 4" for Zero elevation, 3.75" for Anti-Thrust	Height: 4"
18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb
E2607, E2608	E2624, 2625	E2603, E2604



Comfort M2 cushion

May need coccyx relief
 May need postural alignment
 May need support to the greater trochanters



Comfort M2 cushion with ATI

May need coccyx relief
 May need postural alignment
 May need support to the greater trochanters
 May need dynamic positioning



Vicair® Adjuster O2 cushion (standard)

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
 Have normal or impaired sensation
 Experience discomfort from prolonged sitting
 May require different positioning of the left and right femur

Ideal for individuals who

Dimensions

Weight limit

Medicare codes

Width: 18" – 24"	Width: 18" – 24"	Width: 14" – 24"
Depth: 14" – 24"	Depth: 14" – 24"	Depth: 14" – 22"
Height: 4.25" for Zero elevation, 3.75" for Anti-Thrust	Height: 4.25" for Zero elevation, 3.75" for Anti-Thrust	Height: 4"
18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	550 lb
E2607, E2608	E2624, E2625	E2622, E2623



Vicair Adjuster O2 cushion (low)

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
 Have normal or impaired sensation
 Experience discomfort from prolonged sitting
 May require different positioning of the left and right femur

Ideal for individuals who

Dimensions	Width: 16" – 20" Depth: 16" – 20" Height: 2.4" (Kwik Strap only)	Width: 14" – 24" Depth: 14" – 22" Height: 2.4"	Width: 14" – 24" Depth: 14" – 22" Height: 4"
Weight limit	550 lb	440 lb	550 lb
Medicare codes	E2622, E2623	E2624, E2625	E2624, E2625



Vicair Active O2 cushion (low)

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
 Have normal or impaired sensation
 May need postural alignment
 May need support of the greater trochanters
 May require alignment/correction of the pelvis
 May need dynamic positioning
 Those who want easy transfers



Vicair Vector O2 cushion (standard)

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
 Have normal or impaired sensation
 May need postural alignment
 May need support of the greater trochanters
 May require alignment/correction of the pelvis
 May need dynamic positioning



Comfort ACTA-BACK LTS back support

May require a basic amount of postural support
 May require additional lateral trunk support
 May require laterals that are width and height adjustable

Ideal for individuals who

Dimensions	Width: 20" – 24" Length: 14" – 22"	Width: 20" – 24" Length: 14" – 22"	Width: 20" Length: 10" – 22"
Weight limit	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb
Medicare codes	E2620, E2621	E2613, E2614	E2615



Comfort ACTA-EMBRACE® back support

Complex user that truly needs a custom fit



Comfort ACTA-RELIEF® back support

May require a significant amount of postural support
 Those who have moderate to severe kyphosis
 May need the maximum amount of trunk support



Comfort ACTA-BACK® back support

May require a basic amount of postural support
 Those who are more independent with transfers
 Those who have good trunk control

Ideal for individuals who

Dimensions	Width: 20" – 24" Length: 14" – 22"	Width: 20" – 24" Length: 14" – 22"	Width: 20" – 24" Length: 16" – 22"
Weight limit	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb
Medicare codes	E2613, E2614	E2615, E2616	E2620, E2621




Comfort ACTA-BACK Deep back support

May require a significant amount of postural support
 Those who are more dependent with transfers
 May need the maximum amount of trunk support



Comfort ACTA-BACK Contour back support

May require lower thoracic trunk support



Not seeing your size?

The Comfort Inception cushion (E2609) and Custom back (E2617) can be built up to 28" wide. The Inception cushion can also be made thicker as needed.

Contact your Permobil rep for more details.

©2022 Permobil. The following are trademarks and registered trademarks of Permobil: ROHO®, Hybrid Elite®, CONTOUR SELECT®, QUADTRO SELECT®, HIGH PROFILE®, MID PROFILE®, ENHANCER®, ACTA-BACK®, ACTA-EMBRACE®. Vicair® is a registered trademark of Vicair.

Permobil is the distributor of Vicair O2 cushions.

Permobil has a policy of continual product improvement and reserves the right to amend this document. The current version of this document is available at permobil.com.

Medicare Coding: Decisions about how to code a product in the U.S. are the responsibility of the provider/supplier. Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.



Terry Giles

*Custom Fabrication Technician—Power Products
U.S. Navy—Machinery Repairman First Class
(Surface Warfare Specialist/Seabee Combat Warfare Specialist)*

Proudly built in America

Behind every Permobil product is the hard work of hundreds of dedicated Americans. From the world's most advanced power and manual wheelchairs, to proven seating and positioning solutions, each Permobil product in this brochure is assembled with passion and care by employees

in one of our U.S. production locations. Permobil is proud to employ many of our nation's finest veterans and active military personnel who each day continue our mission that EVERYONE has the right to freedom, independence, and mobility.