Power wheelchairs M-Series


T





THE REAL PROPERTY.

Battery status and travel range estimates help you know how far you can go



0

ul 🗢 🛙

Last 7 days

 \otimes

permobil

MyPermobil

MyPermobil connectivity is an available option on all new M-Series models

Available on the App Store

Google Play

Integrated map illustrates daily and historic travel behavio

More adventure

The new Permobil battery and VoltPro charger work together to provide you with up to 25% greater driving distance, reduced charging times, and longer battery life. Use the MyPermobil app to get even more from every adventure with seat activity tracking, battery range estimates, and an integrated map.

4

Reduced

Driving range up to 20 miles

charging times





Longer battery life

App Store is a trademark of Apple Inc. Google Play is a trademark of Google LLC.

More freedom

Corpus®, our most capable power wheelchair seating system, now provides you with even more positioning options. Enhanced stability on the new M-Series enables you to optimize your repositioning, functional reach, and power drive positions for even greater freedom of movement.

Up to 3.2 mph elevated drive speed

Full suspension at any elevation, seat position, or drive speed



Reach for more

ActiveHeight[™] and ActiveReach[™] combine seamlessly to give you the ultimate advantage in accessing your environment. Seat elevation with integrated forward tilt creates the most innovative, functional, and stable seating solution on the market today.

| 12" | 14" |
|-----|------------------|
| | |
| 20° | 20° |
| - | 20° ActiveHei |

3.5° rearward shift during elevate for added stability and reduced footprint

3 points of contact for added stability

Conquer more

FlexLink[™] and FlexLink Pro[™] suspension balances stability, traction, and ride comfort for optimal performance and advanced vibration reduction. Factory optimized and adjustable suspension adapts to your needs and environment. Be ready to handle any adventure with high-powered LED lights to guide your way.

0

0

 \mathbf{O}

FlexLink Pro

FlexLink M3 Corpus

> 2X more powerful LED lights for true pathway lighting and enhanced visibility



M3 Corpus



M5 Corpus

Take control with high-performance motors enabling a top speed of 7.5 mph and exceptional low-end torque that lets you pivot like a pro. Experience unsurpassed stability and ride comfort with FlexLink Pro independent drive wheel suspension. The M5 Corpus will take you to new heights with increased stability and full suspension even in 14" of seat elevation.



BodiLink on Corpus

permobil

Specially designed for the Corpus seating system, BodiLink accessories feature best-in-class hardware. Offering true anatomical support with the highest level of adjustability and comfort to meet your individual positioning needs. Head Support— Easy to adjust piano hinge and double-taper joints for maximum strength

Accent options





Midnight Patriot

American Liberty





Carbon Fiber Mossy Oak[®] Breakup[®] Country™





Water Lily

Fireball



CDG Poly

Please note: Actual product colors may vary from colors shown



Lateral Pelvic/Thigh Support— Easily position along the UniTrack rail with 360° pad adjustment and userfriendly quick release lever

1

-17

0.0

REFERE

| Midnight | |
|----------|--|
| Arctic | |
| Violet | |
| Sky | |

Aqua

Spring

Canary

Sunset

Chili

Flamingo

Platinum

Bronze

Gold



Specifications

| | M3 Corpus | M5 Corpus |
|--------------------------------|---|---|
| Maximum user weight | 330 lb | 330 lb |
| Maximum speed | 6 mph | 7.5 mph |
| Range | 16 mi (Grp 34) 20 mi (Grp 24)* | 20 mi* |
| Base width | 24" | 25" |
| Base length | 35.5" | 37.5" |
| Base turning radius | 20.5" | 21.5" |
| Ground clearance | 3" | 3" |
| Recommended battery type | Group 34 or 24 | Group 24 |
| Weight including batteries | 385 lb (Grp 34) 405 lb (Grp 24) | 423 lb |
| Drive electronics | R-net 120A | R-net 120A |
| Seat-to-floor height | 17.5", 18.5", 19.5" | 17.5", 18.5", 19.5" |
| Seat elevation | Optional 12" | Optional 14" |
| Drive wheel | 14" | 14" |
| Casters | 7" | 8" |
| Tilt options | Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° (<i>up to 300 lb</i>) | Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° <i>(up to 300 lb)</i> |
| Legrest angle - power | 85° – 180° | 85° – 180° |
| Backrest angle - power recline | 85° – 180° | 85° – 180° |
| Backrest height | 20", 23" – 28" (1" incr.) | 20", 23" – 28" (1" incr.) |
| Seat depth | 14" – 22" (1" incr.) | 14" – 22" (1" incr.) |
| Seat width | 17" – 23" (2" incr.) | 17" – 23" (2" incr.) |
| | | |

 * Actual driving range will vary based on driving and battery conditions.

Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions. Chair weights are approximate and include base, batteries, power seat functions, and seating system without accessories.

permobil